

BM3VET Touch Operation Manual

Veterinary Patient Monitor

For Veterinary Use Only



Ver. 1.3

2016.12.22

Table of Contents

Table of Contents	2
1. BASIC	6
1.1 CE Standard Information	7
1.2 Read before Use	9
How to Contact Us	9
Warranty Period	10
Warning, Caution, Note	10
General Precaution on Environment	11
General Precaution on Electric Safety	16
Equipment Connection.....	16
Maintenance and Washing Equipment Connections.....	23
1.3 Product Components	25
Product Outline	25
Principal features of Product.....	25
Product Configuration.....	26
Option Product.....	26
Features of Main Body.....	27
1.4 Function and Key.....	33
External Function.....	33
Operation Key.....	33
1.5 Standard Power Supply Application.....	35
1.6 Battery Power Supply Application.....	36
1.7 General Menu Operation	40
Screen Composition	40
Menu Selection	42
Menu Composition.....	43
2. ANIMAL/DATA MANAGEMENT.....	47
2.1 ADMIT.....	48
ANIMAL TYPE	49
CHANGE ANIMAL INFO.....	49
DEFAULT SETTING.....	50
HEIGHT.....	50
WEIGHT.....	51
2.2 ALARM	52
Alarm for the Product	53
PARAMETER1 LEVEL.....	54
PARAMETER2 LEVEL.....	54
ARRHYTHMIA LEVEL.....	55
ALARM REVIEW.....	55
ALARM VOLUME	57
ALARM PRINT	58
NURSE CALL.....	58
3. SETUP	61

BM3VET Touch Operation Manual

3.1 SETUP	62
SET PARA.....	63
UNIT SELECT	63
USER SERVICE.....	64
SYSTEM.....	71
NETWORK	72
MAKER SERVICE	75
FREEZE MENU	75
4. TREND	76
4.1 TREND	77
GRAPHIC TREND.....	78
TABULAR TREND	79
TREND WINDOW SETUP.....	80
5. ECG	82
5.1 Introduction	83
Colors and Standards of Cables	83
Position of ECG Connector and Measuring Cable.....	83
Attaching Electrodes to the Animal	86
Choosing an ECG lead for Arrhythmia Monitoring.....	87
Information on the ECG waveform.....	87
5 Position of 5-Lead.....	88
Position of 3-Lead Wire Electrodes	89
5.2 ECG Data Window	90
5.3 ECG Data Setup	92
ALARM LIMIT	92
LEAD SELECT	93
QRS VOLUME.....	93
DISPLAY.....	94
ARRHYTHMIA SETTING	97
ST/PVC.....	98
6. SpO₂	104
6.1 Outline	105
SpO ₂ Connector Location and Measuring Cable.....	105
6.2 SpO₂ Data Window	106
Signal and Data Validity	107
6.3 SpO₂ Data Setup	108
ALARM	109
RATE VOLUME	109
LEAD FAULT Condition	110
SPO ₂ Messages.....	110
7. RESPIRATION	111
7.1 Outline	112
7.2 Respiration Data Window	114
7.3 Respiration Data Setup	115
ALARM	116
APNEA DETECT	116
RESPIRATION SPEED	117

BM3VET Touch Operation Manual

RESPIRATION SIZE	117
8. NIBP	118
8.1 Outline.....	119
8.2 NIBP Data Window	121
8.3 NIBP Data Setup.....	123
ALARM	124
CUFF SIZE	125
INFLATION	125
INTERVAL	126
NIBP STAT	127
NIBP VITAL SIGN	128
UNIT SELECT	129
9. EtCO2.....	130
9.1 Introduction.....	131
9.2 EtCO2 Parameter Window	137
9.3 EtCO2 Parameter Setting Menu	138
ALARM LIMIT	139
WAVEFORM SCALE.....	140
EtCO2 SWEEP SPEED	140
APNEA DETECT	141
MODULE INFO	142
MODULE SETUP	143
ZERO.....	145
MODULE RESET	146
9.4 TROUBLESHOOTING	147
10. TEMPERATURE	149
10.1 Outline.....	150
10.2 Temperature Data Window	151
10.3 Temperature Data Setup.....	152
ALARM	152
UNIT SELECT	153
11. PRINT.....	154
11.1 Print.....	155
Printer and Heat Sensitivity Paper	155
Function and Setup Menu	156
11.2 Paper Change	162
12. MESSAGE LIST	163
13. DEFAULT SETTING VALUE	164
13.1 HORSE-ICU Mode	164
13.2 DOG-ICU Mode	167
13.3 PUPPY-ICU Mode	169

BM3VET Touch Operation Manual

13.4 CAT-ICU Mode.....	171
14. TROUBLE SHOOTING.....	173
14.1 Noise in ECG	173
14.2 SpO2 malfunction	174
14.3 Temp malfunction.....	174
14.4 NIBP malfunction	175
14.5 Abnormality in NIBP measurements.....	175
14.6 Failure in battery recharge	176
14.7 Power failure	177
14.8 Periodic noises	178
14.9 Print failure	179
15. SPECIFICATION.....	180
Ease of use	181
Additional Function	181
Monitor Environmental Specifications	181
Power	181
Specification	182
Accessories Included	184
Option	184
Abbreviations and Symbols.....	185
PRODUCT WARRANTY.....	188

1. BASIC

1.1 CE Standard Information

1.2 Read before Use

- Warranty Period
- Warning, Caution, Note
- General Precaution on Environment
- General Precaution on Electric Safety
- Equipment Connection
- Maintenance and Washing Equipment Connection

1.3 Product Components

- Product Outline
- Principal features of Product
- Product Configuration
- Optional Products
- Features of Main Body

1.4 Functions and Key

- External Function
- Operation Key

1.5 Standard Power Supply Application

1.6 Battery Power Supply Application

1.7 General Menu Operation

- Screen Composition
- Menu Selection
- Menu Composition

1.1 CE Standard Information

Electromechanical safety standards met:

Information supplied by the manufacturer of medical devices

1. **EN 60601-1(2006)**
Medical electrical equipment Part1: General requirements for safety
2. **EN 60601-1-2 (2007) (IEC 60601-1-2)**
Electromagnetic Compatibility Requirement and tests
3. **EN 55011:2007+A2:2007 Group 1 Class B(CISPR11) (EN 55011:2009/A1:2010)**
Limits and methods of measurement of radio disturbance characteristics of industrial, scientific and medical (ISM) radio-frequency equipment
4. **IEC 60601-1-4:1996+A1:1999 (EN 60601-1:2006)**
Part 1-4 General requirements for safety Collateral standard: Programmable electrical medical system
5. **IEC 60601-1-6:2010**
Part 1-6 General requirements for safety Collateral standard: Usability
6. **IEC 60601-1-8:2006 (EN 60601-1-8:2007)**
Part 1-8 General requirements for safety Collateral standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
7. **EN 60601-2-30:2000 (IEC 80601-2-30:2013)**
Part 2: Particular requirements for the safety, including essential performance of automatic cycling non-invasive blood pressure monitoring equipment
8. **EN 60601-2-49:2001 (IEC 60601-2-49:2001)**
Part 2: Particular requirements for the safety of multifunction patient monitoring equipment
9. **EN 12470-4:2000+A1:2009**
Performance test for Temperature Clinical thermometers – Part 4: Performance of electrical thermometers for continuous measurement
10. **EN 1060-1:1995+A2:2009, EN 1060-3:1995+A2:2009: (EN ISO 81060-1:2012)**
Performance test for NIBP Non-invasive sphygmomanometers- Part 1: General requirements, Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems
11. **EN ISO 14971:2012**
Medical devices - Application of risk management to medical devices
12. **EN ISO 9919:2009 (ISO 80601-2-61:2011)**
Particular requirements for the basic safety and essential performance of pulse oximeter equipments for medical use.

BM3VET Touch Operation Manual

13. **EN ISO 21647:2009 (ISO 80601-2-55:2011)**
Medical electrical equipment -- Particular requirements for the basic safety and essential performance of respiratory gas monitors
14. **EN 980:2008 (EN ISO 15223-1:2012)**
Symbols for use in the labeling of medical devices (Medical devices -- Symbols to be used with medical device labels, labeling and information to be supplied -- Part 1: General requirements)
15. **EN 1041:2008**
Information supplied by the manufacturer of medical devices
16. **IEC 60601-2-27:2006** : ECG Test
Medical electrical equipment - Part 2-27: Particular requirements for the safety including essential performance, of electrocardiographic monitoring equipment
17. **EN ISO 9919:2005**
Medical electrical equipment -- Particular requirements for the basic safety and essential performance of pulse oximeter equipment for medical use
18. **EN 60601-2-34:2000** : IBP test
Medical electrical equipment – Part2: Particular requirements for the safety, including essential performance, of invasive blood pressure monitoring equipment.

BM3VET Touch Operation Manual

1.2 Read before Use

BIONET services are always available to you.

The following are address and phone numbers for contacting information, services, and product supplies.

How to Contact Us

Product Supply Information

Bionet Co.,Ltd.

5F, Shinsegae I&C Digital Center 61 Digital-ro 31 gil, Guro-gu, SEOUL 08375, REPUBLIC OF KOREA

Tel: +82-2-6300-6410

Fax: +82-2-6499-7789

E-mail: Sales@ebionet.com

Service@ebionet.com

URL: <http://www.ebionet.com>

US Distributor

Bionet America, Inc.

2691 Dow Ave. Ste B

Tustin, CA 92780, USA

Toll Fee: 1-877-924-6638

Tel:1- 714-734-1760

Fax: 1-714-734-1761

www.bionetus.com

sales@bionetus.com

support@bionetus.com

※ In the event of malfunction or failure, contact us along with the model name, serial number, and product name of the equipment.

※ If you need the supply circuit diagram, component list, description and calibration instruction etc. you can contact us we will provide you with it.

BM3VET Touch Operation Manual

Warranty Period

- This product is manufactured and passed through strict quality control and through inspection.
- Compensation standard concerning repair, replacement, refund of the product complies with “Consumer’s protection law” noticed by Korea Fair Trade Commission.
- Warranty period is 4 years.(Four years in USA).
- We will repair or replace any part of the BM3VET Touch found to be defective in usual operating circumstance for free to you.
- This warranty does not apply to any defects caused by improper use, misuse or abuse

Warning, Caution, Note

For special emphasis on agreement, terms are defined as listed below in user’s manual. Users should operate the equipment according to all the warnings and cautions. Indicated in this manual In order to improve the product specifications and features are subject to change without notice.

Warning

To inform that it may cause serious injury or death to the patient, property damage, or material losses

Caution






To inform that it may cause no loss in life but lead to injury

Note

To inform that it is not dangerous but important “note” sign for proper installation, operation, and maintenance of the equipment.

General Precaution on Environment

- Do not keep or operate the equipment in the environment listed below.

	<p>Avoid placing in an area exposed to moisture. Do not touch the equipment with wet hands.</p>		<p>Avoid exposure to direct sunlight</p>
	<p>Avoid placing in an area where there is a high variation of temperature. Operating temperature ranges from 10(C to 40(C). Operating humidity ranges from 30% to 85%.</p>		<p>Avoid placing in the vicinity of Electric heater</p>
	<p>Avoid placing in an area where there is an excessive humidity rise or ventilation problem.</p>		<p>Avoid placing in an area where there is an excessive shock or vibration.</p>
	<p>Avoid placing in an area where chemicals are stored or where there is danger of gas leakage.</p>		<p>Avoid inserting dust or especially metal material into the equipment</p>
	<p>Do not disjoint or disassemble the equipment. This will void your warranty..</p>		<p>Power off when the equipment is not fully installed. Otherwise, equipment could be damaged.</p>

CAUTIONS

Before Installation

Compatibility is critical to safe and effective use of this device. Please contact your local sales or service representative prior to installation to verify equipment compatibility.

Defibrillator Precaution


Patient signal inputs labeled with the CF and BF symbols with paddles are protected against damage resulting from defibrillation voltages. To ensure proper defibrillator protection, use only the recommended cables and lead wires.

Proper placement of defibrillator paddles in relation to the electrodes is required to ensure successful defibrillation.

Disposables

Disposable devices are intended for single use only. They should not be reused as performance could degrade or contamination could occur.

Disposal of your old appliance



1. When this crossed out wheeled bin symbol is attached to a product it means the product is covered by the European Directive 2002/96/EC.
2. All electrical and electronic products should be disposed of separately from the municipal waste stream via designated collection facilities appointed by the government or the local authorities.
3. The correct disposal of your old appliance will help prevent potential negative consequences for the environment and human health.
4. For more detailed information about disposal of your old appliance, please contact your city office, waste disposal service or the shop where you purchased the product.

WARNING

This product contains a chemical known to the State of California to cause cancer, birth defects, or other reproductive harm.

BM3VET Touch Operation Manual

Electrocute Precautions

To prevent skin burns, apply electrocute electrodes as far as possible from all other electrodes, a distance of at 15 cm/6 in. is recommended.

EMC

Magnetic and electrical fields are capable of interfering with the proper performance of the device. For this reason make sure that all external devices operated in the vicinity of the monitor comply with the relevant EMC requirements. X-ray equipment or MRI devices are possible sources of interference as they may emit higher levels of electromagnetic radiation. Also, keep cellular phones and other telecommunication equipment away from the monitor.

BM3VET Touch Operation Manual

CAUTIONS

Intended Use

This device is designed to be used for monitoring the biological vital signs of Canine and Feline and horses. Main functions of the product include displaying information such as ECG, respiration, SpO₂, NIBP, carbon dioxide (CO₂) and temperature on its LCD screen and monitoring parameter, and alarming. It also prints out waves and parameters via a printer.

Application Environment

This device is for use by trained veterinary personnel in veterinary centers. The device is restricted to be used on one patient at a time.

Operator Requirement

Only veterinary personnel who have read the Operator's Manual should use this monitor

Instruction for Use

For continued safe use of this equipment, it is necessary that the instructions are followed. However, instructions listed in this manual in no way supersede established medical practices concerning patient care.

Loss of Data

Should the monitor at any time temporarily lose patient data, the potential exists that active monitoring is not being done. Close Animal observation or alternate monitoring devices should be used until monitor function is restored.

If the monitor does not automatically resume operation within 60 seconds, power cycle the monitor using the power on/off switch. Once monitoring is restored, you should verify correct monitoring state and alarm function.

Maintenance

Regular preventive maintenance should be carried out annually (Technical inspections). You are responsible for any requirements specific to your country and locality.

MPSO

The use of a multiple portable socket outlet (MPSO) for a system will result in an enclosure leakage current equal to the sum of all individual earth leakage currents of the system if there is an interruption of the MPSO protective earth conductor. Do not use an additional extension cable with the MPSO as it will increase the chance of the single protective earth conductor interruption.

Negligence

BIONET does not assume responsibility for damage to the equipment caused by improperly vented cabinets, improper or faulty power, or insufficient wall strength to support equipment mounted on such walls.

NOTES

Power Requirements

Before connecting the device to the power line, check that the voltage and frequency ratings of the power line are the same as those indicated on the unit's label. If this is not the case, do not connect the system to the power line until you adjust the unit to match the power source. In U.S.A, if the installation of this equipment will use 240V rather than 120V, the source must be a center-tapped, 240V, single-phase circuit.

Restricted Sale

U.S.A federal law restricts this device to sale by or on the order of a licensed veterinarian..

Supervised Use

This equipment is intended for use under the direct supervision of trained veterinary personnel in veterinary centers. The device is restricted to be used on one patient at a time.

Ventilation Requirements

Set up the device in a location which affords sufficient ventilation. The ventilation openings of the device must not be obstructed. The ambient conditions specified in the technical specifications must be ensured at all times.

·Put the monitor in a location where you can easily see the screen and access the operating controls.

·This product is protected against the effects of cardiac defibrillator discharges to ensure proper recovery, as required by test standards. (the screen may blank during a defibrillator discharge but recovers within second as required by test standards.)

Reference Literature

Medical Device Directive 93/42/EEC

EN 60601-1/1990 +A1: 1993 +A2 : 1995 : Medical electrical equipment.

General requirements for safety

EN 60601-1-1/9. 1994 +A1 12.95: General requirements for safety.

General Precaution on Electric Safety

Warning

Check the items listed below before operating the equipment.

1. Be sure that AC power supply line is appropriate to use. (AC100 - 240V)
2. Be sure that the power source is the one supplied from Bionet. (DC18V,2.8A, BPM050S18F02 Made in BridgePower Co., Ltd.)
3. Be sure that the entire connection cable of the system is properly and firmly fixed.
4. Be sure that the equipment is completely grounded. (If not, there might be problems in the product.)
5. The equipment should not be placed in the vicinity of electric generators, X-ray, broadcasting apparatus to eliminate electrical noise during operation. Otherwise, it may cause incorrect results.

Note

The Equipment should be placed far from generators, X-ray equipment, broadcasting equipment or transmitting wires, so as to prevent electrical noise from being generated during operation, When these devices are near the Equipment, it can produce inaccurate measurements. For BM3VET TOUCH both independent circuit and stable grounding are essentially required. In the event that the same power source is shared with other electronic equipment, it can also produce inaccurate output.

Warning

Do not make contact with the Animal while operating the machine It may cause serious danger to the users. Use only the provided cables.
A warning that other cables and accessories may negatively affect EMC performance

Warning

In case the Equipment does not operate as usual or damaged, do not use on Animal, and contact to the medical equipment technician of the hospital or the equipment supply division.

Note

BM3VET TOUCH is classified as follows:
- BM3VET TOUCH classifies as Class I, BF & CF concerning electric shock. It is not proper to operate this Equipment around combustible anesthetic or dissolvent.
- Noise level is B class regarding IEC/EN 60601-1 and the subject of Nose is B level concerning IEC/EN60601-1-2.

Equipment Connection

BM3VET Touch Operation Manual

For measurements in or near the heart we recommend connecting the monitor to the potential equalization system. Use the green and yellow potential equalization cable and connect it to the pin

labeled with the symbol .

Manufacturer's declaration - electromagnetic emission

The BM3VET TOUCH system is intended for use in the electromagnetic environment specified below. The customer or the user of BM3VET TOUCH system should ensure that it is used in such an environment		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The BM3VET TOUCH system uses RF energy only for its internal function. Therefore. Its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The BM3VET TOUCH system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supplies buildings used for domestic purposes.
Harmonics emission IEC 61000-3-2	A	
Voltage fluctuation IEC 61000-3-3	Complies	

BM3VET Touch Operation Manual

Manufacturer's declaration - electromagnetic immunity

The **BM3VET TOUCH** system is intended for use in the electromagnetic environment specified below. The customer or the user of the **BM3VET TOUCH** system should ensure that it is used in such an environment

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic Environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV Contact 8 kV Air	6 kV Contact 8 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast Transient / burst IEC 61000-4-4	2kV for power supply lines 1kV for input/output lines	2kV for power supply lines 1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	1 kV differential mode 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) Magnetic field IEC 61000-4-8	3.0 A/m	3.0 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short Interruptions and Voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5cycle 40% U_T (60% dip in U_T) for 5 cycle 70% U_T (30% dip in U_T) for 25 cycle <5% U_T (<95% dip in U_T) for 5 s	<5% U_T (>95% dip in U_T) for 0.5cycle 40% U_T (60% dip in U_T) for 5 cycle 70% U_T (30% dip in U_T) for 25 cycle <5% U_T (<95% dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BM3VET TOUCH system requires continued operation during power mains interruptions, it is recommended that the BM3VET TOUCH system be powered from an uninterruptible power supply or a battery

Note: U_T is the a.c. mains voltage prior to application of the test level.

The **BM3VET TOUCH** system is intended for use in the electromagnetic environment specified below

The customer or the user of the **BM3VET TOUCH** system should ensure that it is used in such an environment

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
---------------	----------------------	------------------	---------------------------------------

BM3VET Touch Operation Manual

<p>Conducted RF IEC 61000-4-6</p>	<p>3 Vrms 150 kHz to 80 MHz</p>	<p>3 Vrms 150 kHz to 80 MHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the BM3VET TOUCH system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p style="text-align: center;">Recommended separation distance</p> $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$
<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80.0 MHz to 2.5 GHz</p>	<p>3 V/m 80.0 MHz to 2.5 GHz</p>	<p style="text-align: center;">Recommended separation distance</p> $d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,</p> <p>(a) Should be less than the compliance level in each frequency range (b).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>
<p>Note 1) U_T is the A.C. mains voltage prior to application of the test level.</p>			
<p>Note 2) At 80 MHz and 800 MHz, the higher frequency range applies.</p>			
<p>Note 3) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the EUT is used exceeds the applicable RF compliance level above, the EUT should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the EUT.</p>			
<p>b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than $[V_1]$ V / m.</p>			
<p>Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the BM3VET Touch system.</p>			

BM3VET Touch Operation Manual

The **BM3VET TOUCH** system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the **BM3VET Touch** system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BM3VET TOUCH system as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power (W) of transmitter	Separation distance (m) according to frequency of transmitter		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

BM3VET Touch Operation Manual

Immunity and Compliance Level			
Immunity test	IEC 60601 Test Level	Actual Immunity Level	Compliance Level
Conducted RF IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 Vrms, 150 kHz to 80 MHz	3 Vrms, 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	3 V/m, 80 MHz to 2.5 GHz	3 V/m, 80 MHz to 2.5 GHz	3 V/m, 80 MHz to 2.5 GHz

Guidance and manufacturer's declaration - electromagnetic immunity

The **BM3VET TOUCH** system is intended for use in the electromagnetic environment specified below. The customer or the user of the **BM3VET TOUCH** system should ensure that it is used in such an environment

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms 150 kHz to 80 MHz	BM3VET TOUCH system must be used only in a shielded location with a minimum RF shielding effectiveness and, for each cable that enters the shielded location with a minimum RF shielding effectiveness and, for each cable that enters the shielded location
Radiated RF IEC 61000-4-3	3 V/m 80.0 MHz to 2.5 GHz	3 V/m 80.0 MHz to 2.5 GHz	Field strengths outside the shielded location from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3V/m. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Note 2) It is essential that the actual shielding effectiveness and filter attenuation of the shielded location be verified to assure that they meet the minimum specification.

a- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength outside the shielded location in which the EUT is used exceeds 3V/m, the EUT should be observed to verify normal operation.

If abnormal performance is observed, additional measures may be necessary, such as relocating the EUT or using a shielded location with a higher RF shielding effectiveness and filter attenuation.

Note

For Type A Professional ME Equipment intended for use in domestic establishment instructions for use includes a warning:

This ME equipment is intended for use by professional healthcare personnel only.

BM3VET Touch Operation Manual

Caution

In the hospital, doctors and Animals are exposed to dangerous, uncontrollable compensating currents. These currents are due to the potential differences between connected equipment. The safety solution to the problem is accomplished with EN60601-1;1996.

Biocompatibility

When used as intended, the parts of the product described in this operator manual, including accessories that come in contact with the Animal during the intended use, fulfill the biocompatibility requirements of the applicable standards. If you have questions about this matter, please contact BIONET or its representatives.

BM3VET Touch Operation Manual

Maintenance and Washing Equipment Connections

Various methods can be used to clean the BM3VET TOUCH and its accessories. Please follow the methods mentioned below to avoid unnecessary damage or contamination to the Equipment.

We do not repair free of charge regardless of warranty period if it is contaminated or damaged by using dangerous material not designated for washing.

Cleaning Applied Parts

Do not permit any liquid to enter the monitor case and avoid pouring it on the monitor while cleaning. Do not allow water or cleaning solution to enter the connectors of jack cover.

Recommended cleaning agents:

Alcohol (Ethanol 70%, Isopropanol 70%, Window cleaner)

Ammonias (Dilution of ammonia <3%, Window cleaner)

Tensides (dishwasher detergents) (Edisonite schnellreiniger[®], Alconox[®])

Cables and Leadwires

CAUTION

Do not use acetone or ketone solvents for cleaning; do not use an autoclave or steam cleaner.

Cables and leadwires can be cleaned with a warm, damp cloth and mild soap, or isopropyl alcohol wipes. For more intensive disinfecting (near sterile) Ethylene Oxide (ETO) is acceptable but will reduce the useful lifetime of the cable or leadwire.

CAUTION

The decision to sterilize must be made per your institution's requirements with an awareness of the effect on the integrity of the cable or leadwire.

Note

The Equipment needs a safety inspection once a year. Please refer to user's guide or service manual for the procedure.

Please check carefully both frame and sensor, after cleaning the Equipment, Do not use equipment that is worn out or damaged.

At least once a month, clean and wipe off the frame by using a soft cloth after wetting it with water and alcohol. Do not use lacquer, thinner, ethylene, or oxidizer which may lead to damage to the equipment.

Make sure both cables and accessories are free of dust or contaminants, and wipe them off with soft cloth wetted with warm water (40°), and at least once a week, clean them by using clinical alcohol. Do not submerge the accessories under any liquid or detergent. Also, make sure liquids do not penetrate into the Equipment or probe.

Disinfecting

Do not mix disinfecting solutions (such as bleach and ammonia) as hazardous gases may result.

Clean equipment before disinfecting.

Recommended disinfecting agents:

Aldehyde based (Cidex[®] activated dialdehyde solution, Gigasept)

Alcohol base (Ethanol 70%, Isopropanol 70%, Spitacid[®], Streilium fluid[®], Cutasept[®], Hospisept[®], Tinktur forte, Sagrosept[®], Kodan[®])

BM3VET Touch Operation Manual

Caution

Always follow all local laws and recommendations for disposal of single use and/or contaminated items.

Caution

There is back-up battery inside system. When users dispose of this battery, please follow all local laws and recommendations. .

Warning

Check the electrodes of batteries before changing them.

- Operate BM3VET TOUCH with internal electric power supply when unsure of external ground connection.
- Remove the 1st Battery when not using equipment for an extended period of time to avoid any damage.

For other applied parts such as temperature sensors, pulse oximetry probes, and NBP cuffs, you must consult the manufacturer for cleaning, sterilization, or disinfecting methods.

1.3 Product Components

Product Outline

BM3VET TOUCH monitor is a product used for monitoring biological information of Canine and Feline. Main functions of the product include displaying information such as ECG, respiration, SpO2, NIBP, EtCO2 and temperature on its LCD screen and monitoring parameters, and alarming. It also prints out waves and parameters via a printer.

Principal features of Product

BM3VET TOUCH is a small-size multifunctional monitoring unit for an Animal designed for easy usage during movement. It features a DC power supply (Bridgepower, BPM050S18F02, DC 18V, 2.8A). The equipment also measures major parameters such as ECG, respiration rate, SpO2, pulse rate, NIBP, EtCO2, and temperature, displaying them on an 8-inch color LCD screen. It also enables users to check waves and parameters and other vital signs of an Animal via the 58mm thermal printer and monitor the Animal using the alarm system. With B-Link software, up to 128 hours of saved parameter data can be transferred to a Windows based computer through an Ethernet network.

Warning

Use only the accessories provided by us. Otherwise, Animal and user may be exposed to danger.

Warning

BEFORE USE — Before putting the system into operation visually inspect all connecting cables for signs of damage. Damaged cables and connectors must be replaced immediately. Before using the system, the operator must verify that it is in correct working order and operating condition. Periodically, and whenever the integrity of the product is in doubt, test all functions.

BM3VET Touch Operation Manual

Product Configuration

1. Main body of BM3VET Touch Monitor	1 EA
2. 3-Lead vet ECG Cable (3CBL-400, 3WIRE-430)	1 EA
3. 3-Lead vet Extension Cable	1 EA
4. NIBP extension tube (NBPCBL-400)	1 EA
5. NIBP vet cuff infant reusable	1 EA
6. SpO ₂ sensor extension cable (SPCBL-400)	1 EA
7. Reusable multisite SpO ₂ probe	1 EA
8. DC Adaptor (BPM050S18F02 made in Bridgepower Co., Ltd.)	1 EA
9. Operator`s Manual	1 EA
10. Chart Paper (PAPER-400)	2Roll
11. Temperature probe	1 EA
12. Printer(Built in)	1 EA
13. Esophageal ECG/temperature probe	1 EA
14. Transflectance SpO ₂ sensor	1 EA

Option Product

1. 5-Lead Animal Cable (MECA5(AHA),MECE5(IEC))
2. EtCO₂ Module
3. Li-ion Battery (2150mAh, 10.8V)
4. Sidestream EtCO₂ Module (Respironics)
5. Mainstream EtCO₂ Module (Respironics)
6. Sidestream EtCO₂ airway adapter sampling kit
7. Mainstream EtCO₂ airway adapter

Warning

In order to avoid electrical shock, do not open the cover. Disassembling of the equipment should be done only by service personnel authorized by BIONET


Warning

Users must pay attention on connection of any auxiliary device via LAN port or nurse calling. Always consider about summation of leakage current, please check if the auxiliary device is qualified by IEC 60601-1, or consult your hospital biomedical engineer.


BM3VET Touch Operation Manual

Features of Main Body

FRONT




BACK




BM3VET Touch Operation Manual

Right Side



Left Side



BM3VET Touch Operation Manual

Accessories





	3-Lead ECG Extension Cable
	5-Lead ECG Extension Cable
	3-Lead vet ECG Cable
	5-Lead vet ECG Cable
	Reusable multisite SpO2 Probe
	SPO2 extension cable (2m)
	NIBP Extension Tube (3m)
	Reusable NIBP Infant Cuff Cuff Size : 210 * 60 Range : 8 to 13 cm
	Temperature Probe (Surface/Skin)
	Temperature Probe (Rectal/Esophageal)

BM3VET Touch Operation Manual



Equipment Sign

	<p>ATTENTION : Consult accompanying documents</p>
	<p>TYPE CF APPLIED PART : Insulated (floating) applied part suitable for intentional external and internal application to the Animal including direct cardiac application. "Paddles" outside the box indicate the applied part is defibrillator proof. Medical Standard Definition : F-type applied part(floating/insulated) complying with the specified requirements of IEC 60601-1/UL 2601-1/CSA 601.1 Medical Standards to provide a higher degree of protection against electric shock than that provided by type CF applied parts.</p>
	<p>TYPE BF APPLIED PART : Insulated (floating) applied part suitable for intentional external and internal application to the Animal excluding direct cardiac application. "Paddles" outside the box indicate the applied part is defibrillator proof. Medical Standard Definition : F-type applied part (floating/insulated) complying with the specified requirements of IEC 60601-1/UL 2601-1/CSA 601.1 Medical Standards to provide a higher degree of protection against electric shock than that provided by type BF applied parts.</p>

BM3VET Touch Operation Manual

	Ground
	Printer
	Serial Port
	LAN Port
	AUX Connector Port
	DC Input Indicator
	DSUB 15pin external VGA port
	Battery Operation Indicator
	WIRELESS LAN power output Port
	DC Input Connector

BM3VET Touch Operation Manual


	TOUCH SCREEN LOCK
	NIBP
T	Temperature
F	Function
	Power on
	Power off
	Respiration
	ECG
	Heart Pulse

BM3VET Touch Operation Manual

1.4 Function and Key

External Function

The front panel of this product consists of an LCD screen and five function keys and one trim knob.



Operation Key

1. Power : Switches on and off the Power.
2. Function Key – Alternates between display modes.
3. Blood Pressure : Manually completes measuring blood pressure.
4. Printer : Prints out the waves selected from the menu until the key is pressed to stop.
5. Alarm : Stop alarm sound.
 - First press stops the current alarm for one minute
 - Second press stops the all alarm for five minutes.
 - Third press will stops the all alarms.
 - Forth press resets alarms back to the original setting.
6. Trim Knob : This key is used to select menu by turning it clock or anticlockwise to move cursors.
7. Alarm + Function: Touchscreen, key, rotary wheel lock function on and off.

BM3VET Touch Operation Manual

Lock : Press the alarm and function key at the same time until lock icon is displayed.


Unlock : Press the alarm and function key at the same time until Unlock icon is displayed.




lock icon



unlock icon



Use knob or touch screen to select menu



BM3VET Touch Operation Manual

1.5 Standard Power Supply Application

DC Power

- Product information
- Manufacture: Bridgpower corp.
- Model name: BPM050S18F02
- Input power:

Rated Voltage 100 – 240V

Rated Line Frequency 50 – 60 Hz

Current 1.5A Max at 100 VAC input

Protection Internal primary current Fuse (2.0A)


- Output power:

Capable of supplying +18VDC at 2.8A

Voltage +18VDC +/- 5%

Over current Protection 3.36 – 5.6A

DC Power LED is lighted on when the DC Power is plugged into the inlet on the back of the product. A press of power key makes the machine ready for use.



Warning

This equipment must only be connected to a supply mains with ground. Noise or distortion of signals using non-off-the-shelf products rather than adapters supplied by our company may be caused.

BM3VET Touch Operation Manual

1.6 Battery Power Supply Application


Battery power can be supplied for enabling portable use or for use during DC power failure.

Operation

1. Battery Power LED is lighted on when the machine is in use.
2. The DC/battery power is only sustainable for 1 and a half hours.
3. Battery is automatically charged when the machine is connected to DC Power Supply. Battery LED is lighted on after blinking.
4. The charging status of the batteries is displayed with 5 green boxes, e.g. charging. (0% -> 25% -> 50% -> 75% -> 100%)


- Battery: 031PpTC(3ICR19/65)(10.8V, 2150mAh/23.22Wh)

The Lithium-Ion battery is a rechargeable battery containing Lithium-Ion cells. Each battery contains an integrated electronic fuel gauge and a safety protection circuit.



BM3VET Touch Operation Manual

5. The discharge condition of battery is indicated with on of 5 yellow boxes, each box showing a different level of charge available.
(100% -> 75% -> 50% -> 25%)



When remaining battery is less than 25%, the battery icon box is turned to blinking red. The device will be turned off automatically after 5 minutes from that warning sign. In case of that warning sign with blinking red icon, charge the device immediately with DC power adaptor which is provided from BIONET.




- Battery charging time: More than 6 hours
- Continuous battery use time: Lowest 1 hour to highest 2 hours continuous use (buffering)

Warning

Check the electrodes of batteries before charging them.

6. Battery status indication: When battery is disconnected from equipment or out of order, it is shown by a red 'X' as shown below.



7. Low power supply: When power is less than 16V, the battery indication disappears and the "LOW" indication is active.




Display of low power

Note

When the batteries are replaced, remove and replace the DC adapter.

To insert and remove the battery pack.

Assembly or replacement, as shown in the figure below.



The Impact of Lithium-Ion Battery Technology on the Battery

The following are the key points you should know about Lithium-Ion battery technology:

The battery will discharge on its own, even when it is not installed in a monitor. This discharge is the result of the Lithium-Ion cells and the bias current required for the integrated electronics.

By the nature of Lithium-Ion cells, the battery will self-discharge.

The self-discharge rate doubles for every 10°C (18°F) rise in temperature.

The capacity loss of the battery degrades significantly at higher temperatures.

As the battery ages, the full-charge capacity of the battery will degrade and be permanently lost. As a result, the amount of charge that is stored and available for use is reduced.

BM3VET Touch Operation Manual

Conditioning Guideline

The battery in the monitor should be fully charged and discharged every six months and condition it using the battery charger.

Storage Guideline

Store the battery outside of the monitor at a temperature between 20°C to 25°C (68°F to 77°F).

When the battery is stored inside a monitor that is powered by an AC power source, the battery cell temperature increases by 15°C to 20°C (59°F to 68°F) above the room's ambient temperature. This reduces the life of the battery.

When the battery is stored inside a monitor that is continuously powered by an AC power source and is not powered by battery on a regular basis, the life of the battery may be less than 12 months.

BIONET recommends that you remove the battery and store it near the monitor until it is needed for transport..

How to Recycle the Battery

When the battery no longer holds a charge, it should be replaced. The battery is recyclable. Remove the old battery from the monitor and follow your local recycling guidelines.

WARNING


EXPLOSION HAZARD —

DO NOT incinerate the battery or store at high temperatures. Serious injury or death could result.

BM3VET Touch Operation Manual

1.7 General Menu Operation

Screen Composition



- Real Time Wave Window : Displays measured results by up to seven waves.
- Menu icon windows: A menu displays icons for each function.
- Menu Select Window : Menus appear when they are activated..
- Parameter Window : Measured and setup data are displayed in five windows.


The screen consists of a total of two modes

NORMAL MODE: As the figure shows the waveform parameters of the screen


PARAMETER MODE: Selected figures show only the 4 parameters of the screen

BM3VET Touch Operation Manual

NORMAL MODE




PARAMETER MODE



BM3VET Touch Operation Manual

Menu Selection




When the Trim Knob Key is turned, menus are selected in the order indicated above. The above screen shows that the MORE menu is selected. The menu move to the right in the order of MORE MENU → ECG → NIBP → SpO₂ → RESP (EtCO₂) → TEMP. An inactivated window is jumped off.

Menu Composition


More Menu Window

When the additional menu is selected it will set and cancel the functions.



Numerical value sign widow


This window displays a measured parameter, function setup, and the boundary of parameter values.



BM3VET Touch Operation Manual


Menu selection by using Trim Knob key

As the key is turned to the right, the menu selection moves clockwise. As the key is turned to the left, the menu selection moves counterclockwise. The menu selection is activated when you depress Trim Knob key.



Menu selection with touch keys

Touch the desired menu, the menu will be selected.



BM3VET Touch Operation Manual

Word feature menu

The following figure shows the screen where the word sequence menu is activated within the word sequence correction menu. Here, the cursor moves over the words when the Trim Knob key is turned in the clockwise direction.

Touch the desired menu box, select the menu is available



Word feature menu


The following figure shows the screen where the word sequence menu is activated within the word sequence correction menu. Here, the cursor moves over the words when the Trim Knob key is turned in the clockwise direction. To enter letters and numbers at the touch of their letters and numbers after the 'SET' button is pressed



BM3VET Touch Operation Manual

Operation menu

The setup value changes without a selection when the menu is moved.



2. ANIMAL/DATA MANAGEMENT

2.1 ADMIT

ANIMAL TYPE
CHANGE ADMIT INFO
DISCHARGE
HEIGHT
WEIGHT

2.2 ALARM

Alarm for the Product
PARAMETER1 LEVEL
PARAMETER2 LEVEL
ARRHYTHMIA LEVEL
ALARM REVIEW
ALARM VOLUME
ALARM PRINT
NURSE CALL

2.1 ADMIT
CHANGE ANIMAL INFO
ANIMAL TYPE
DEFAULT SETTING
HEIGHT UNIT
WEIGHT UNIT

The screenshot shows a vertical menu titled "ADMIT" with a close button (X) in the top right corner. The menu contains five options, each in a separate box:

- ANIMAL TYPE: HORSE
- CHANGE ANIMAL INFO
- DEFAULT SETTING
- WEIGHT UNIT: KG
- HEIGHT UNIT: CM

The bottom portion of the menu is a solid grey area.

BM3VET Touch Operation Manual

ANIMAL TYPE

Set the patient environment of equipment in Animal Type.

HORSE : LARGE ANIMAL // DOG: MEDIUM ANIMAL

PUPPY : SMALL ANIMAL // CAT: TINY ANIMAL

ADMIT

ANIMAL TYPE: HORSE

CHANGE ANIMAL INFO

DEFAULT SETTING

WEIGHT UNIT: KG

HEIGHT UNIT: CM

HORSE

DOG

PUPPY

CAT

CHANGE ANIMAL INFO

Last and first name (11 letters for each), sex (male or female), date of birth, weight, height, and Animal ID (11 characters)

ADMIT

ANIMAL TYPE: HORSE

CHANGE ANIMAL INFO

DEFAULT SETTING

WEIGHT UNIT: KG

HEIGHT UNIT: CM

LAST NAME: JOHN

FIRST NAME: WASHINGTON

ANIMAL ID: 0012198367752

SEX: MALE

BIRTH DATE: 03-06-1981

AGE: 31


HEIGHT: 178.0 Cm

WEIGHT: 80.0 Kg

BM3VET Touch Operation Manual


DEFAULT SETTING

Parameter range settings, alarm settings, and Animal-specific initialization settings.



HEIGHT

Unit of height is set as Cm / Inch.



BM3VET Touch Operation Manual

WEIGHT

Unit of weight is set as Kg / LBS.

ADMIT X

ANIMAL TYPE: HORSE

CHANGE ANIMAL INFO

DEFAULT SETTING

WEIGHT UNIT: KG

HEIGHT UNIT: CM

This screenshot shows a vertical menu with a blue header bar containing the text 'ADMIT' and a close button 'X'. Below the header are five menu items, each in a grey box with a black border. The 'WEIGHT UNIT: KG' item is highlighted with a blue border.

ADMIT X

ANIMAL TYPE: HORSE

CHANGE ANIMAL INFO

DEFAULT SETTING

WEIGHT UNIT: LBS

HEIGHT UNIT: CM

This screenshot is identical to the one on the left, but the 'WEIGHT UNIT: LBS' item is highlighted with a blue border.

2.2 ALARM

Alarm is divided into two, alarm for the Animal's condition and for the product's condition. The Animal's alarm sounds when the diagnostic functions (ASYSTOLE, VTAC/VFIB, and VTAC) are detected. Each alarm sound differs in order and volume according to the levels of HIGH, MEDIUM, LOW and MESSAGE.

HIGH	 -5			
MEDIUM	 -3			
LOW	 -1			
MESSAGE				



: Alarm sounds



: Number flashes



: Waves are printed out

Alarm for the Product

The machine gives alarm sounds for its system with a related message flashing.

LOW



-1



ALARM LIMITS : The machine enables one to see and change the limits of alarm for all parameter functions.


ALARM PRINT : with an ON/OFF setup, the related information is printed out whenever an alarm is given.

ALARM VOLUME : volume of each alarm can be adjusted in 10 steps.

ALARM LEVEL : Priority of each parameter alarm can be set up.

ALARM REVIEW : Shows the priority order information for all alarms of each measurement.

NURSE CALL : Set the feature of the NURSE CALL.



Able to see and change all the alarm functions.

BM3VET Touch Operation Manual

PARAMETER1 LEVEL

ECG, NIBP, SpO2, RESP and TEMP information about all of the alarm settings.

PARAMETER1 LEVEL	HR MESSAGE 150 50	NIBP-S MESSAGE 150 50	SPO2-% MESSAGE 150 50	RESP MESSAGE 150 50
PARAMETER2 LEVEL	ST MESSAGE 150 50	NIBP-M MESSAGE 150 50	SPO2-PR MESSAGE 150 50	RESP-A MESSAGE 150 50
ARRHYTHMIA LEVEL	PVC MESSAGE 20 0	NIBP-D MESSAGE 150 50	NIBP-PR MESSAGE 160 60	TEMP1 MESSAGE 42.0 30.0
ALARM REVIEW	LEAD FAULT MESSAGE	LOW BATTERY MESSAGE		

ALARM LEVEL	ALARM
MESSAGE	ON
MESSAGE	HIGH
LOW	50
MEDIUM	LOW
HIGH	150

1	2	3	-
4	5	6	CLR
7	8	9	SET
0	.	<-	

PARAMETER2 LEVEL

EtCO2 information about all of the alarm settings.

PARAMETER1 LEVEL	EtCO2 MESSAGE 150 50	FiCO2 MESSAGE 150 50	AWRR MESSAGE 30 10	APNEA LOW 20 0
PARAMETER2 LEVEL				
ARRHYTHMIA LEVEL				
ALARM REVIEW				


ALARM LEVEL	ALARM
MESSAGE	OFF
MESSAGE	HIGH
LOW	50
MEDIUM	LOW
HIGH	150

1	2	3	-
4	5	6	CLR
7	8	9	SET
0	.	<-	

BM3VET Touch Operation Manual


ARRHYTHMIA LEVEL

Diagnostics when the alarm is set.



ALARM REVIEW

After an alarm is triggered the alarms and data wave pattern can be reviewed. Set up the priority of each parameter alarm.



BM3VET Touch Operation Manual

ALARM LIST

When an alarm activates, this shows the order of the alarms. Up to 20 cases of alarm can be stored.

ALARM [X]

PARAMETER1 LEVEL	ALARM LIST	LIST	TIME	KIND	UP
		ASYSTOLE	2011/03/18 10:22:53	MEDIUM	[Vertical Scroll Bar]
		VTAC	2011/03/18 10:22:53	HIGH	
		VTAC/FIB	2011/03/18 10:22:53	HIGH	
		BIGEMINY	2011/03/18 10:22:53	HIGH	
		ACC VENT	2011/03/18 10:22:53	HIGH	
		SPO2	2011/03/18 10:22:53	MEDIUM	
		RESP	2011/03/18 10:22:53	MESSAGE	
		IRRGULAR	2011/03/18 10:22:53	HIGH	
		IRRGULAR	2011/03/18 10:22:53	HIGH	
		ACC VENT	2011/03/18 10:22:53	HIGH	
		PVC	2011/03/18 10:22:53	MEDIUM	DN

PARAMETER1 LEVEL
 PARAMETER2 LEVEL
 ARRHYTHMIA LEVEL
 ALARM REVIEW
 ALARM VOLUME: OFF
 ALARM PRINT: ON
 NURSE CALL:

SAVE CONDITION: HIGH
 MESSAGE
 LOW
 MEDIUM
 HIGH

ALARM [X]

PARAMETER1 LEVEL	RETURN	2011/03/24 07:22:10	MEDIUM
	80BPM	[ECG Waveform]	
	98% SpO2	[SpO2 Waveform]	
	20RPM	[Respiration Waveform]	
	NIBP : 120/80/(91)	[NIBP Waveform]	
	T1 : 36.5	[Temperature Waveform]	
	IBPT : 120/81/(92)	[IBPT Waveform]	
	22:21:10	22:21:13	22:21:17

PARAMETER1 LEVEL
 PARAMETER2 LEVEL
 ARRHYTHMIA LEVEL
 ALARM REVIEW
 ALARM VOLUME: OFF
 ALARM PRINT: ON
 NURSE CALL :

[Horizontal Scroll Bar]

BM3VET Touch Operation Manual

SAVE CONDITION

This determines the alarm level of parameters which are saved in the alarm list when alarm occurs. If a higher level of alarm occurs than the previously determined alarm level, data would be saved in the alarm list.

The screenshot shows the 'ALARM' window with a sidebar on the left containing settings for PARAMETER1 LEVEL, PARAMETER2 LEVEL, ARRHYTHMIA LEVEL, ALARM REVIEW, ALARM VOLUME: OFF, ALARM PRINT: ON, and NURSE CALL:. The main area is divided into 'ALARM LIST' and 'LIST'. The 'ALARM LIST' section has a 'SAVE CONDITION : HIGH' button highlighted in blue, along with 'MESSAGE', 'LOW', 'MEDIUM', and 'HIGH' buttons. The 'LIST' section is a table with columns for LIST, TIME, and KIND. The table contains the following data:

LIST	TIME	KIND
ASYSTOLE	2011/03/18 10:22:53	MEDIUM
VTAC	2011/03/18 10:22:53	HIGH
VTAC/FIB	2011/03/18 10:22:53	HIGH
BIGEMINY	2011/03/18 10:22:53	HIGH
ACC VENT	2011/03/18 10:22:53	HIGH
SPO2	2011/03/18 10:22:53	MEDIUM
RESP	2011/03/18 10:22:53	MESSAGE
IRRGULAR	2011/03/18 10:22:53	HIGH
IRRGULAR	2011/03/18 10:22:53	HIGH
ACC VENT	2011/03/18 10:22:53	HIGH
PVC	2011/03/18 10:22:53	MEDIUM

ALARM VOLUME

Set the alarm volume at 10 levels.

The screenshot shows the 'ALARM' window with the 'ALARM VOLUME: OFF' setting highlighted in blue. The main area displays a grid of buttons for setting the alarm volume at various levels:

OFF	10%	20%
30%	40%	50%
60%	70%	80%
90%	100%	

BM3VET Touch Operation Manual

ALARM PRINT

ON / OFF settings for when the alarm information is printed on thermal paper.

The screenshot shows the 'ALARM' settings window. On the left side, there is a vertical list of settings: PARAMETER1 LEVEL, PARAMETER2 LEVEL, ARRHYTHMIA LEVEL, ALARM REVIEW, ALARM VOLUME: OFF, ALARM PRINT: ON (highlighted with a blue border), and NURSE CALL:. The main area of the window is a large grey rectangle.

NURSE CALL

When an alarm is triggered, this activates the NURSE CALL function.

The screenshot shows the 'ALARM' settings window with the 'NURSE CALL:' option highlighted. The settings are as follows: PARAMETER1 LEVEL, PARAMETER2 LEVEL, ARRHYTHMIA LEVEL, ALARM REVIEW, ALARM VOLUME: OFF, ALARM PRINT: ON, and NURSE CALL: (highlighted). The NURSE CALL: option is expanded to show sub-options: NURSE CALL TYPE: NORMAL OPEN, NURSE CALL DURATION: CYCLING (with sub-options ONE TIME and CONTINUE), and NURSE CALL LEVEL: LOW (with sub-options LOW, MEDIUM, and HIGH).


BM3VET Touch Operation Manual

NURSE CALL TYPE

NURSE CALL function call when an alarm condition is met.

NORMAL OPEN: RELAY OPEN when ALARM does not ring, CLOSE when ALARM does ring.

NORMAL CLOSE: RELAY CLOSE when ALARM does not ring, OPEN when ALARM does ring.




NURSE CALL DURATION
NURSE CALL alarm calls when the

situation is set to output mode.

ONE TIME: After ALARM occurs, set the RELAY to be ON for 3 seconds then OFF

CYCLING: Relay will cycle between ON and OFF in 1-second intervals.

CONTINUE: After ALARM occurs, set the RELAY to be ON for 60 seconds then OFF.



BM3VET Touch Operation Manual

NURSE CALL LEVEL

NURSE CALL alarm level is set to operate.

LOW : If the alarm is raised above the LOW level NURSE CALL call.

MEDIUM: If the alarm is raised above the MEDIUM level NURSE CALL call.

HIGH: If the alarm is raised above the HIGH level NURSE CALL call.

The screenshot shows the 'ALARM' settings screen. On the left is a vertical menu with the following options: PARAMETER1 LEVEL, PARAMETER2 LEVEL, ARRHYTHMIA LEVEL, ALARM REVIEW, ALARM VOLUME: OFF, ALARM PRINT: ON, and NURSE CALL:. The main area contains the following settings:

- NURSE CALL: ON
- NURSE CALL TYPE: NORMAL OPEN
- NURSE CALL DURATION: CYCLING
- ONE TIME (button)
- CONTINUE (button)
- CYCLING (button)
- NURSE CALL LEVEL: LOW (highlighted with a blue border)
- LOW (button)
- MEDIUM (button)
- HIGH (button)

3. SETUP

3.1 SETUP

SET PARA
UNIT SELECT
USER SERVICE
SYSTEM
NETWORK
MAKER SERVICE
FREEZE MENU

3.1 SETUP



The Settings menu is displayed after pressing the icon shown above.


SET PARA : On and Off measurement parameters of patient monitor.

UNIT SELECT: Select a unit for pressure, ST, and temperature.

USER SERVICE: This is the menu to set the connection used to interface with an external computer.

SYSTEM: Display system version information.


MAKER SERVICE: This is the basic adjustment menu used to adjust the features of this product.



BM3VET Touch Operation Manual

SET PARA

Select measurement function to use.



UNIT SELECT


This is the menu for converting the units of BM3VET TOUCH.

The units of parameters for pressure, ST LEVEL, Temperature are able to convert

Pressure: kPa \leftrightarrow mmHg

ST mm \leftrightarrow mV


Temperature: $^{\circ}\text{C}$ \leftrightarrow $^{\circ}\text{F}$



BM3VET Touch Operation Manual

USER SERVICE


The user is able to set the communication parameters, power supply filter.



BM3VET Touch Operation Manual

SET UNIT NAME

Set up for Equipment name




BM3VET Touch Operation Manual

SET BED NUMBER

Set up for Animal bed number.

Allowable settings are from 1 to 10.




BM3VET Touch Operation Manual

SET DATE & TIME

Set date and time of equipment.


Press the SET button after each input change you want to change the year, month, day, hour, minute, and second item during the setting will be entered.



BM3VET Touch Operation Manual


W-LAN

Power supplying of W-LAN module could be adjusted with this function



KEY SOUND


This is the menu for KEY SOUND to ON/OFF.



BM3VET Touch Operation Manual


DEMO

Set ON/OFF DEMONTRATION of equipment.



AC FILTER

AC FILTER is function where you can set power supply frequency. This feature is required because power supply frequency can be different from one country to another. . (The selectable frequencies are 50Hz and 60Hz, OFF.)




BM3VET Touch Operation Manual

COLOR SELECT

This is the menu to set the waveform and parameter color selection.

It has ten colors below table.


The color of parameter could be changed to one of ten colors from following table.



BM3VET Touch Operation Manual

SYSTEM

System able to change and verify Equipment version information and system information.



VGA OUTPUT: Turn on VGA output in order to connect external monitor.

BM3VET Touch Operation Manual

NETWORK

Setup network connection.

The screenshot shows the 'SETTING' window with a sidebar menu on the left containing 'SET PARA', 'UNIT SELECT', 'USER SERVICE', 'SYSTEM', 'NETWORK', and 'MAKER SERVICE'. The 'NETWORK' option is selected. The main area contains several settings:

CENTRAL: ON	DHCP: ON
HOST IP	192.168.030.100
DEVICE IP	192.168.030.101
SUBNET	255.255.255.000
GATEWAY	192.168.030.001
MAC ADDR	00:02:BD:80:00:00

CENTRAL : If this is ON, system transfers parameter data to the computer where Host IP is assigned.

DHCP: If this is ON, system automatically obtains the device's IP from router.

HOST IP: a computer's IP address where B-Link software was installed.

DEVICE IP: Patient monitor's IP address (This will be assigned automatically, if DHCP is ON).

SUBNET: This will be assigned automatically, if DHCP is ON(Default is 255.255.255.0).

GATEWAY : This will be assigned automatically, if DHCP is ON (For more GATEWAY information, contact your network administrator).

The screenshot shows the 'HOST IP' input screen. At the top, the text 'HOST IP' is displayed. Below it, the IP address '192.168.30.101' is shown in a sequence of input fields, with '192' currently selected and highlighted. Below the input fields is a numeric keypad with buttons for digits 1-9, 0, a decimal point, and a left arrow. There are also 'CLR' and 'SET' buttons.

BM3VET Touch Operation Manual

DEVICE IP: Press the SET button to set the address of the sending and receiving equipment.

DEVICE IP

192 . 168 . 30 . 101

1 2 3 CLR
4 5 6
7 8 9 SET
0 . ←

SUBNET : Press the SET button to set.

SUBNET

255 . 255 . 255 . 000

1 2 3 CLR
4 5 6
7 8 9 SET
0 . ←

BM3VET Touch Operation Manual

GATEWAY : Press the SET button to set the address that set up a connection at the network settings window




MAC ADDR : Press the SET button to set the hardware address for the network settings.
(MAC ADDR ≡ Impossible to modify by the user.)



MAKER SERVICE

Maker service is a menu used by the manufacturer.



FREEZE MENU

If you select the icon which is located in the far left of the icon menu by controlling the rotary switch, the wave window is held and is maintained as the previous status, at the same time the parameter windows is normally showing the current Animal's status.

Whenever selecting the FREEZE menu, the FREEZE and RELEASE are repeated by turns.



The FREEZE is released by the following two conditions.

1. 3 minutes after selecting FREEZE menu.
2. Selection of the release FREEZE menu

Note

Unlike regular screen freeze screen prints during waveform output parameters of the state and its parameters.

4. TREND

4.1 TREND

GRAPHIC TREND
TABLE TREND
TREND WINDOW SETUP

BM3VET Touch Operation Manual

4.1 TREND

TREND shows saved data graphically displayed with numeric values.

Real-time data recording interval is 1 minute. Amount of saving time for this data will be 168hours. The following entries are stored.

(HR,ST,PVC,NIBP(S/M/D),RESP,SPO2%,SPO2-PR,TEMP, EtCO2,FICO2,AWRR)

TREND						
TABULAR TREND	GRAPHIC TREND	TREND WINDOW SETUP				
TABULAR TREND						
	15-FEB 2011					22:13
	15-FEB 22:12	15-FEB 22:12	15-FEB 22:12	15-FEB 22:12	15-FEB 22:12	15-FEB 22:12
HR	80	80	80	80	80	80
ST	0.2	0.2	0.2	0.2	0.2	0.2
PVC	0	0	0	0	0	0
NIBP-S	120	120	120	120	120	120
NIBP-D	80	80	80	80	80	80
NIBP-M	93	93	93	93	93	93
NIBPPR						
RESP						
SPO2-%	99	99	99	99	99	99
SPO2-R	80	80	80	80	80	80
TEMP1	36.5	36.5	36.5	36.5	36.5	36.5
AWRR	20	20	20	20	20	20
EtCO2	32	32	32	32	32	32
FICO2	0	0	0	0	0	0
ALARM						


X : Move to main screen.

0.5 1 1.5 3 6 : Time period set menu
1 5 15 30 60

BM3VET Touch Operation Manual

GRAPHIC TREND

Wave Data can be stored and reviewed according to section.



TIME PERIOD

One can set up and store data and time that one can see on a screen.



BM3VET Touch Operation Manual

TABULAR TREND

One can see the stored data at the time previously set up.

TABULAR TREND		15-FEB 2011 22:13				
	15-FEB 22:12	15-FEB 22:12	15-FEB 22:12	15-FEB 22:12	15-FEB 22:12	15-FEB 22:12
HR	80	80	80	80	80	80
ST	0.2	0.2	0.2	0.2	0.2	0.2
PVC	0	0	0	0	0	0
NIBP-S	120	120	120	120	120	120
NIBP-D	80	80	80	80	80	80
NIBP-M	93	93	93	93	93	93
NIBP-PR						
RESP	20	20	20	20	20	20
SPO2-%	99	99	99	99	99	99
SPO2-R	80	80	80	80	80	80
TEMP1	36.5	36.5	36.5	36.5	36.5	36.5
AWRR	20	20	20	20	20	20
EtCO2	32	32	32	32	32	32
FiCO2	0	0	0	0	0	0
ALARM						

TIME INTERVAL



One can store data and set up time.



BM3VET Touch Operation Manual

TREND WINDOW SETUP


Set the trend display window that will show on the real time wave window.



BM3VET Touch Operation Manual


SET TREND PARA

Set visible parameters on the screen.



TIME PERIOD PARA

Set visible time period in a screen.



TREND PRINT

Graphic: select the number which selects a graphic trend and press print to print the selected trend.
Table: select the table number to be print and press print to print all the data in the selected Animal admit (Admit) table.

5. ECG

5.1 Introduction

- Colors and Standards of Cables
- Position of ECG Connector and Measuring Cable
- Attaching Electrodes to the Animal
- Choosing an ECG lead for Arrhythmia Monitoring
- Information on the ECG waveform
 - 5 Position of 5-Lead
- Position of 3-Lead Wire Electrodes

5.2 ECG Data Window

5.3 ECG Data Setup

- TRACE 1 LEAD SELECT
- ALARM LIMIT
- ALARM
- QRS VOLUME
- ECG SIZE
- HEART RATE SOURCE
- ECG SPEED
- ANALYSIS SETTING

BM3VET Touch Operation Manual

5.1 Introduction

It calculates the heart rate with 3 or 5 leads ECG signal acquisition and performs the alarms according to the set values.

Colors and Standards of Cables

AHA : American Heart Association (U.S.A. Certification)

IEC : International Electro technical Commission (Europe Certification)

3LEAD / 5LEAD


Leadwire	AHA Color code	AHA Label	IEC Color code	IEC Label
Right arm	White	RA	Red	R
Left arm	Black	LA	Yellow	L
Right leg	Green	RL	Black	N
Left leg	Red	LL	Green	F
V1(precordial)	Brown	V1	White	C1

AHA : American Heart Association (U.S.A. standard)

IEC : International Electro technical Commission (Europe standard)

Position of ECG Connector and Measuring Cable

ECG connector +detect cable



BM3VET Touch Operation Manual



IEC 3LEAD CABLE



AHA 3LEAD CABLE



IEC 5LEAD CABLE



AHA 5LEAD CABLE

BM3VET Touch Operation Manual



IEC 3LEAD



AHA 3LEAD



IEC 5LEAD



AHA 5LEAD

Attaching Electrodes to the Animal

1. Shave excess hair. With a piece of cotton pad moistened with alcohol, clean the Animal's skin where the electrodes should be mounted. Avoid wrinkled or uneven skin areas. Wipe off the alcohol with a dry cotton pad.
2. Open the electrode package and take out the electrode.
3. Remove the backing paper from the electrode. Be careful not to touch the adhesive side.
4. Attach the disposable electrode to the previously cleaned skin. Avoid wrinkled and uneven skin areas.
5. The electrode lead is connected to the monitor onto the electrode.
6. Fasten the electrode lead to the skin with surgical tape with an extra length of wire between the tape and the electrode. This prevents body movement from moving the electrode lead.

Note

- ✓ To maintain good contact between the electrode and skin, check that the paste of the disposable electrode is not dry.
- ✓ When contact of the disposable electrode becomes poor, replace the electrode with a new one immediately. Otherwise, contact impedance between the skin and electrode increase and the correct ECG cannot be obtained.
- ✓ If the contact is bad before the expiration date on the package, replace the electrode with a new one.
- ✓ To obtain a stable ECG wave form rub the skin with "skin Pure" skin preparation gel or tincture of Benzion.
- ✓ Use only the CE certified disposable electrode.

Choosing an ECG lead for Arrhythmia Monitoring

It is very important to select a suitable lead for arrhythmia monitoring.
Guidelines for non-paced Animals:


- ✓ QRS should be tall and narrow(recommended amplitude $> 0.5\text{mV}$)
- ✓ R wave should be above or below the baseline (but not bi-phasic)
- ✓ T wave should be smaller than $1/3$ R-wave height.
- ✓ The P-wave should be smaller than $1/5$ R-wave height.

For paced Animals, in addition to the above,:

- ✓ Not wider than the normal QRS
- ✓ The QRS complexes should be at least twice the height of pace pulses.
- ✓ Large enough to be detected, with no re-polarization.


To prevent detection of P-waves or baseline noises as QRS complexes, the minimum detection level for QRS complexes is set at 0.15mV . Adjusting the ECG wave size on the monitor display(gain adjustment)does not affect the ECG signal which is used for arrhythmia analysis. If the ECG signal is too small, you may get false alarms for asystole.

Information on the ECG waveform




When ECG signal is 80bpm T-wave duration is 180ms, and the QT interval is 350ms.


5 Position of 5-Lead



Position of 3-Lead Wire Electrodes



5.2 ECG Data Window



Note

ECG Wave Display is always on when the cable is connected.

The heart rate is calculated by a moving average. The monitor detects 8 consecutive beats, averages the R-R intervals of the latest 8 beats and uses this average to calculate the current heart rate. When a new beat is detected, the heart rate is recalculated using the latest 8 beats. The heart rate display is updated every 3 seconds.

Heart rate meter updates a new heart rate for a step increase or decrease in 10 seconds maximum.

When ventricular tachycardia is detected, the alarm set in 5 seconds maximum.

Check that the delay time of the output signal (alarm trigger 80ms maximum) is within the range of the connected equipment.

BM3VET Touch Operation Manual

Safety Precautions

Warning

CABLES — Route all cables away from Animal's throat to avoid possible strangulation.

CONDUCTIVE CONNECTIONS — Extreme care must be exercised when applying medical electrical equipment. Many parts of the human/machine circuit are conductive, such as the Animal, connectors, electrodes, transducers. It is very important that these conductive parts do not come into contact with other grounded, conductive parts when connected to the isolated Animal input of the device. Such contact would bridge the Animal's isolation and cancel the protection provided by the isolated input. In particular, there must be no contact of the neutral electrode and ground.

DEFIBRILLATION — Do not come into contact with Animals during defibrillation. Otherwise serious injury or death could result.

To avoid the risk of serious electrical burn, shock, or other injury during defibrillation, all persons must keep clear of the bed and must not touch the Animal or any equipment connected to the Animal.

After defibrillation, the screen display recovers within 10seconds if the correct electrodes are used and applied in accordance with the manufacturer's instructions.

Animal cables can be damaged when connected to a Animal during defibrillation. Check cables for functionality before using them again.

The peak of the synchronized defibrillator discharge should be delivered within 60ms of the peak of the R wave. The signal at the ECG output on the Animal monitors is delayed by a maximum of 30ms.

If the ECG waveform on the screen is too unstable to synchronize with the Animal's heart beat because of the following reason, remove the cause of an alarm, message, or unstable ECG, and then use a stable ECG lead for synchronization.

- ✓ ECG electrode is detached or broken. Lead wire is detached or broken.
- ✓ Lead wire moves. AC interference, EMG noise or noise from ESU is superimposed.
- ✓ Connection cable is broken or has a short circuit. Connector has poor contact.

INTERFACING OTHER EQUIPMENT — Devices may only be interconnected with each other or to parts of the system when it has been determined by qualified biomedical engineering personnel that there is no danger to the Animal, the operator, or the environment as a result. In those instances where there is any element of doubt concerning the safety of connected devices, the user must contact the manufacturers concerned (or other informed experts) for proper use. In all cases, safe and proper operation should be verified with the applicable Manufacturer's instructions for use, and system standards IEC 60601-1-1/EN 60601-1-1 must be complied with.


Electrosurgery Unit

- ✓ Electrosurgical units(ESU) emit a lot of RF interference. If the monitor is used with an ESU,RF interference may affect the monitor operation.
- ✓ Locate the monitor as far as possible from the ESU. Locate them on opposite sides of the operating table, if possible.
- ✓ Connect the monitor and ESU to different AC outlets located as far as possible from each other.
- ✓ When using this monitor with an electrosurgical unit, its return plate and the electrodes for monitoring must be firmly attached to the Animal. If the return plate is not attached correctly,it may burn the Animal's skin where the electrodes are attached.

BM3VET Touch Operation Manual

5.3 ECG Data Setup


A setup window appears at lower part of the screen when the Trim Knob Key is pressed in the ECG Parameter Window. Selection is made by pressing the Trim Knob Key, while movement across the menu is performed by turning the key either clock or anticlockwise.



ALARM LIMIT

Alarm Limit is 0 ~ 300BPM.

ECG alarm feature ON / OFF and the menu is set to LEVEL.




BM3VET Touch Operation Manual

LEAD SELECT

Select channels from I to V in ECG

Lead I, II, III show up in case of connecting 3-Leads Animal Cable.


Lead I, II, III, aVR, aVL, aVF, V show up in case of connecting 5-Leads Animal Cable.



QRS VOLUME

Move the Key to select a volume rate from OFF, 10% to 100%.


When QRS volume "ON", SpO2 volume setting is set OFF automatically.



BM3VET Touch Operation Manual

DISPLAY

Set the sweep speed and waveform size.



ECG FILTER


One may select from three frequency types for WAVE FILTER.

MONITOR 0.5Hz ~ 40Hz

MODERATE 0.5Hz ~ 25Hz

MAXIMUM 5Hz ~ 25Hz

DIAGNOSIS 0.05Hz ~ 150Hz




BM3VET Touch Operation Manual

ECG SWEEP SPEED


ECG speed on the LCD is 25 mm/s.

Speed is changeable to 6.25, 12.5, 25, 50mm/s.



ECG SIZE

The size is changeable to X0.25, X0.5, X1, X2, X4.




BM3VET Touch Operation Manual

ECG VIEW

The number of ECG waves can be configured with this function.


There are 2 traces of 1 CH data at the ECG wave.

- 3LEAD or 5LEAD : 1CH




HR Source

ECG or SpO2 can be selected as heart rate source.



ARRHYTHMIA SETTING

Analysis setting is divided to 3 menus.



ARRHYTHMIA : Sets up ON/OFF to indicate detection of diagnosis (Asys, VTAC/VFIB and VTAC).

OFF: Do not perform arrhythmia diagnosis.

LETHAL: Performs the detection of Asys, VTAC/VFIB, and VTAC at the selected lead

The Analysis algorithm uses a selected lead I, II, III, or V lead for ECG and arrhythmia analysis.

ASYSTOLE

Ventricular asystole occurs whenever the displayed heart rate drops to zero.

PVC

Isolated premature ventricular complexes occur when a premature ventricular beat is detected and has non-ventricular beats before and after.


VFIB/VTAC

Ventricular fibrillation occurs when the ECG waveform indicates a chaotic ventricular arrhythm.

BM3VET Touch Operation Manual


ARRHYTHMIA ALARM LEVEL

Diagnostic alarm level is set.




ST/PVC

ST signal and setting related ST menu and PVC menu.




BM3VET Touch Operation Manual

ST ANALYSIS: ON/OFF ST analysis signal.




PVC ANALYSIS: Decision maker to display PVC value sign with ON/OFF




BM3VET Touch Operation Manual

INITIAL SETUP: ST measurements to factory settings (ISO R-: 80, ST R+: 108)



ST MEASUREMENT CONDITION fine-tune the ISO and ST in order to position the cursor keys to select the rotary and then to be adjusted and controlled at ISO and ST TOUCH TOUCH button arrow and then fine-tuning is possible when TOUCH.


MEASUREMENT CONDITION: ST measurement condition setting



BM3VET Touch Operation Manual

PACE : Sets up ON/OFF to indicate that the Animal has PACE.

The PACE menu option enables/disables the pacemaker detection program.



Be aware of the warning below when monitor an animal with a pacemaker.

Warning

FALSE CALLS—False low heart rate indicators or false asystole calls may result with certain pacemakers because of electrical overshoots.

MONITORING PACEMAKER ANIMALS—Monitoring of pacemaker Animals can only occur with the pace program activated.

PACEMAKER SPIKE—An artificial pacemaker spike is displayed in place of the actual pacemaker spike. All pacemaker spikes appear uniform. Do not diagnostically interpret pacemaker spike size and shape.

ANIMAL HAZARD—A pacemaker pulse can be counted as a QRS during asystole in either pace mode. Keep pacemaker Animals under close observation.

PACEMAKER ANIMALS. Rate meters may continue to count the pacemaker rate during occurrences of cardiac arrest or some arrhythmias. Do not rely entirely upon rate meter ALARMS. Keep pacemaker Animals under close surveillance.

CAUTION

FDA POSTMARKET SAFETY ALERT

The United States FDA Center for Device and Radiological Health issued a safety bulletin October 14, 1998. this bulletin states "that minute ventilation rate-adaptive implantable pacemakers can occasionally interact with certain cardiac monitoring and diagnostic programmed rate."

The FDA further recommends precautions to take into consideration for Animals with these types of pacemakers. These precaution for Animals with these types of pacemakers. These precautions include disabling the rate responsive mode and enabling an alternate pace mode. For more information contact:

Office of Surveillance and Biometrics, CDRH, FDA
1350 Packard Drive, Mail Stop HFZ-510 Rockville, MD 20850 U.S.A

NOTE

ECG monitoring with Animals in non-invasive trans coetaneous pacemakers may not be possible due to large amounts of energy produced by these devices. Monitoring ECG with an external device may be needed.

WARNINGS

VENTRICULAR ARRHYTHMIAS

The arrhythmia analysis program is intended to detect ventricular arrhythmias. It is not designed to detect a trial or supra ventricular arrhythmias. Occasionally it may incorrect identify the presence or absence of an arrhythmia. Therefore, a physician must analyze the arrhythmia information in conjunction with other clinical findings.

SUSPENDED ANALYSIS

Certain conditions suspend arrhythmia analysis. When suspended, arrhythmia conditions are not detected and alarms associated with arrhythmias do not occur. The messages which alert you to the conditions causing suspended arrhythmia analysis are : ARR OFF, ARRHYSUSPEND, LEADS FAIL, ALARM PAUSE, ALL ALARMS OFF, and DISCHARGED.

Trouble shooting

Problem :

Inaccurate heart rate and/or false a systole.

Solution :

Check ECG signal from Animal:

1. Check/adjust lead placement.
2. Check/perform skin preparation.
3. Check/replace electrodes.

Check amplitude of ECG waveform:

1. Select ECG parameter label.
2. Select DISPLAY LEAD,
3. Scroll through all ECG leads and check for 0.5mV amplitude at normal (1X) size. (at least 0.5mV amplitude is required for QRS detection.) for borderline signals, validate on a graph.
4. If amplitudes are low, electrodes may need to be repositioned or replaced.

Problem :

False ventricular calls.

Solution :

Check ECG signal from Animal: (the chest lead may exhibit polarity changes which may occasionally cause an inaccurate call.)

1. Check/adjust lead placement.
2. Check/perform skin preparation.
3. Check/replace electrodes. (if chest lead is a problem, move the chest lead to another chest position or leg position.)

Problem :

Inaccurate pacemaker detection

Solution :

Use pacemaker processing:

1. Select ECG parameter label.
2. Display the lead of ECG with the greatest amplitude in the top waveform position.
3. Select ANALYSIS SETTINGS.
4. SELECT DETECT PACE.

6. SpO₂

6.1 Outline

SpO₂ Connector Location and Measuring Cable

6.2 SpO₂ Data Window

Signal and Data Validity

6.3 SpO₂ Data Setup

ALARM

RATE VOLUME

LEAD FAULT Condition


SPO₂ Messages

6.1 Outline

SPO2 monitoring is a noninvasive technique used to measure the amount of oxygenated hemoglobin and pulse rate by measuring the absorption of selected wavelengths of light. The light generated in the probe passes through the tissue and is converted into an electrical signal by the photodetector in the probe. The monitor processes the electrical signal and displays on the screen a waveform and digital values for SpO2 and pulse rate. It detects SpO2 in the way of transmitting the red and infrared rays into the capillary vessel to take the pulsation. Also performs the alarm function according to the set values.

SpO2 Connector Location and Measuring Cable


SpO₂ connector



Note

The signal input is a high-insulation port and it is defibrillator proof ()
 The insulated input ensures Animal safety and protects the device du  s fibrillation and electro-surgery.

6.2 SpO2 Data Window



The current SPO₂ value and the derived pulse rate (RATE) are displayed. The block sets indicate the strength of the signal (twenty block bars indicate the strongest signal). The SPO₂ measurements are averaged over a 6-second period of time. The monitor display is updated every second. The SPO₂ monitoring features are found in the SPO₂ menu. These features include alarm limit adjustment, display of RATE, and RATE volume.

Note

SpO₂ WAVE SIZE is changed automatically.

Signal and Data Validity

It is extremely important to determine that the probe is attached to the Animal correctly and the data is verifiable. To make this determination, three indications from the monitor are of assistance—signal strength bar, quality of the SPO₂ waveform, and the stability of the SPO₂ values. It is critical to observe all three indications simultaneously when ascertaining signal and data validity.

Signal Strength Bar

The signal strength bar is displayed within the SPO₂ values window. This bar consists of 20 blocks set depending on the strength of the signal. Proper environmental conditions and probe attachment will help to ensure a strong signal.

Quality of SPO₂ Waveform

Under normal conditions, the SPO₂ waveform corresponds to (but is not proportional to) the arterial pressure waveform. The typical SPO₂ waveform indicates not only a good waveform, but helps the user find a probe placement with the least noise spikes present. The figure below represents an SPO₂ waveform of good quality.



Good Quality SPO₂ Waveform

If noise (artifact) is seen on the waveform because of poor probe placement, the photodetector may not be flush with the tissue. Check that the probe is secured and the tissue sample is not too thick. Pulse rate is determined from the SPO₂ waveform which can be disrupted by a cough or other hemodynamic pressure disturbances. Motion at the probe site is indicated by noise spikes in the normal waveform. (See the figure below.)



SPO₂ Waveform with Artifact

Stability of SPO₂ Values

The stability of the displayed SPO₂ values can also be used as an indication of signal validity. Although stability is a relative term, with a small amount of practice one can get a good feeling for changes that are artifactual or physiological and the speed of each. Messages are provided in the SPO₂ values window to aid you in successful SPO₂ monitoring.

WARNING

In the monitoring of Animals the coincidence of adverse conditions may lead to a disturbed signal going unnoticed. In this situation artifacts are capable of simulating a plausible parameter reading, so that the monitor fails to sound an alarm. In order to ensure reliable Animal monitoring, the proper application of the probe and the signal quality must be checked at regular intervals.


6.3 SpO2 Data Setup

ALARM : SpO₂ alarm limit set up.

RATE VOLUME : SpO₂ volume set up.

HR SOURCE: Heart rate in ECG window source selection menu(same as ECG menu).


SWEEP SPEED: ECG and SpO₂ waveform display sweep speed selection menu(same as ECG menu).



BM3VET Touch Operation Manual


ALARM

Two menus: ALARM LIMIT, ALARM provided in the alarm menu
Number setting of alarm value of %SpO₂ is 0 ~ 100
Warning sound or message displays when an alarm is triggered.



RATE VOLUME

Move the KEY to select the volume from OFF to 100%.
The SpO₂ volume setting turns on a tone which sounds each time an SpO₂ pulse is detected. This is a variable pitch tone which changes as the Animal's saturation level changes.



LEAD FAULT Condition

When using a reusable finger probe, there is a system alarm to alert you when the probe is off the Monitor. The monitor defaults this " LEAD FAULT" condition as a System Warning alarm, however, you can set it as a System ALARM LEVEL in Monitor Defaults.

SPO2 Messages

Below is a list of system status alarm messages which may be displayed in the SPO2 parameter window during monitoring.

CHECK PROBE

Reusable finger probe is off the Animal. Check the probe. *The factory default for this alarm is MESSAGE ALARM.*

PULSE SEARCH

Detection by the monitor of a repeatable pulse has ceased. Check the Animal and the probe site.

POOR SIGNAL

The SPO2 signal is too low. No SPO2 data is displayed. This can be due to a low Animal pulse, Animal motion, or some other interference. Check the Animal and the probe.

LOST SIGNAL

SPO2 data continues to be displayed, but the quality of the signal is questionable. Check the Animal and the probe.

ARTIFACT

It indicates that something happened to the pulses; determine if the artifact to be abnormal and irregular

7. RESPIRATION

7.1 Outline

7.2 RESPIRATION Data Window

7.3 RESPIRATION Data Setup

Alarm
APNEA DETECT
RESPIRATION SPEED
RESPIRATION SIZE


BM3VET Touch Operation Manual

7.1 Outline

Respiration expands the skin area of the chest, causing changes in the resistance of the skin. Through this it calculates respiration value per minutes and performs the alarm function according to the set value.

Respiration Connector and Measuring Cable

Respiration Connector



Respiration Measuring Cable



IEC 3LEAD CABLE



AHA 3LEAD CABLE



IEC 5LEAD CABLE



AHA 5LEAD CABLE

BM3VET Touch Operation Manual



IEC 3LEAD



AHA 3LEAD



IEC 5LEAD




AHA 5LEAD

Note

RR uses the ECG cable and lead wires.

Alcohol (Ethanol 70%, Isopropanol 70%, Window cleaner)
Alcohol (Ethanol 70%, Isopropanol 70%, Window cleaner)

7.2 Respiration Data Window




7.3 Respiration Data Setup

ALARM: Respiration alarm setting menu

APNEA DETECT: A menu to setup APNEA alarm display

SWEEP SPEED: A menu to setup Wave Display of speed


RESP SIZE: A menu to setup Wave Display



BM3VET Touch Operation Manual

ALARM

Alarm menu provide ALARM LIMIT and ALARM SOUND .




Alarm Limit of Respiration Numeric Value is 5 ~ 150bpm

Alarm Limit of RESPIRATION APNEA Numeric Value is 3 ~ 30sec.

Warning sound or message displays when Respiration ALARM occurs.

APNEA DETECT


Deciding function of activating Apnea Alarm



BM3VET Touch Operation Manual


RESPIRATION SPEED

Wave pattern speed is 6.25 , 12.5 , 25 mm/s.



RESPIRATION SIZE

Set wave pattern size X2~ X10.



8. NIBP

8.1 Outline

8.2 NIBP Data Window

8.3 NIBP Data Setup

- ALARM
- CUFF SIZE
- INFLATION
- INTERVAL
- NIBP STAT
- NIBP VITAL SIGN
- UNIT SELECT


BM3VET Touch Operation Manual

8.1 Outline

This function is to measure minimum, Maximum and average blood pressure by using Oscillometric method

Position of NIBP Connector and cuff

NIBP Connector




INFANT CUFF



BM3VET Touch Operation Manual

POSITION OF CUFF (CAT)



POSITION OF CUFF (DOG)



WARNING

Noninvasive blood pressure monitoring is not recommended for Animals with hypotension, hypertension, arrhythmias or extremely high or low heart rate. The software algorithm cannot accurately compute NIBP or Animals with these conditions.

Note

Tubes between the cuff and the monitor should not be kinked or blocked.

The air pad should be exactly over the branchial artery. Tubing is immediately to the right or left of the branchial artery to prevent kinking when elbow is bent.


BM3VET Touch Operation Manual

The maintenance is performed every 2 years.


Check the following list to ensure device operates properly and safety at all times.

1. Check for proper cuff size.
2. Check for residual air left in the cuff from a previous measurement.
3. Make sure cuff is not too tight or too loose.
4. Make sure cuff and heart are at same level, otherwise hydrostatic pressure will offset the NIBP value.
5. Minimize Animal movement during measurement.
6. Watch for pulses paradox us.
7. Check for leak in cuff or tubing.
8. Animal may have a weak pulse.

8.2 NIBP Data Window



BM3VET Touch Operation Manual



POWER OFF

When power is cut off during pressure, air runs out of the CUFF automatically.

8.3 NIBP Data Setup

ALARM : A menu to set the Alarm

CUFF SIZE : A menu to select cuff size


INFLATION: Initial Pressurization setting menu

INTERVAL : A menu to set Interval time when measuring the blood pressure automatically

NIBP STAT : 5 Minutes continuous measurement

NIBP VITAL SIGN : History display of NIBP measurement value

UNIT SELECT: A menu to select the pressure unit




BM3VET Touch Operation Manual

ALARM

The alarm provides ALARM LIMIT and ALARM SOUND.

Alarm setting Numeric Value of Systolic, Diastolic, and mean pressure is 10 ~ 360mmHg.


The menu which decide activate of warning sign and message display when the respiration alarm is on.



BM3VET Touch Operation Manual

CUFF SIZE

The user can select a CUFF between ADULT and NEONATAL.




INFLATION

It is a function for set the maximum initial inflation pressure value. The range of initial inflation pressure value of BM3VET TOUCH is as follows.

LARGE : 120 – 250 mmHg / MEDI. : 120 – 250mmHg / SMALL : 60 – 140mmHg

Default Value: LARGE : 170 mmHg / MEDI. : 140mmHg / SMALL : 120mmHg





BM3VET Touch Operation Manual

INTERVAL

This menu is used for selecting intervals when measuring the blood pressure automatically.

Select a target interval from 1min, 2, 3, 4, 5, 10, 15, 20, 30, 1hour, 2, 4, 8.

INTERVAL is set after the start, press the NIBP KEY to initiate readings.



BM3VET Touch Operation Manual

If you select the icon  from the main screen of the equipment you can select the measurement cycle.




Warning

Periodically check Animal limb circulation distal to the cuff. Check frequently when using auto NBP in 1 and 2 minute intervals. Intervals below 10 minutes are not recommended for extended periods of time.

NIBP STAT

5 minutes to continuous measurement mode.



BM3VET Touch Operation Manual

NIBP VITAL SIGN


15 recently measured blood pressure values and pulse rates are recorded.

NIBP			
	TIME	SYS / DIA (MEAN)	PR
ALARM	2011/03/18 10:22:53	120 / 80 (94)	65BPM
CUFF SIZE: LARGE	2011/03/18 10:23:53	120 / 80 (94)	65BPM
INFLATION: 170mmHg	2011/03/18 10:24:33	120 / 80 (94)	65BPM
INTERVAL: OFF	2011/03/18 10:25:23	120 / 80 (94)	65BPM
NIBP STAT: OFF	2011/03/18 10:26:43	120 / 80 (94)	65BPM
NIBP STAT: OFF	2011/03/18 10:27:53	120 / 80 (94)	65BPM
NIBP STAT: OFF	2011/03/18 10:28:58	120 / 80 (94)	65BPM
NIBP STAT: OFF	2011/03/18 10:29:25	120 / 80 (94)	65BPM
NIBP STAT: OFF	2011/03/18 10:30:28	120 / 80 (94)	65BPM
NIBP VITAL SIGN	2011/03/18 10:31:12	120 / 80 (94)	65BPM
NIBP VITAL SIGN	2011/03/18 10:32:28	120 / 80 (94)	65BPM
UNIT SELECT: mmHg	2011/03/18 10:33:34	120 / 80 (94)	65BPM
UNIT SELECT: mmHg	2011/03/18 10:34:43	120 / 80 (94)	65BPM
UNIT SELECT: mmHg	2011/03/18 10:35:28	120 / 80 (94)	65BPM

BM3VET Touch Operation Manual

UNIT SELECT

It is a function to set blood pressure measurement unit.
The blood pressure measurement unit provides mmHg and kPa.



Warning

Pay attention not to block connecting hose when you put cuff on Animal.
Check cuff or hose connection for leaks periodically. Measurements can be inaccurate if air leaks.

9. EtCO₂

9.1 INTRODUCTION

9.2 EtCO₂ Parameter Window

9.3 EtCO₂ Parameter Setting Menu

ALARM LIMIT

WAVEFORM SCALE

EtCO₂ SWEEP SPEED

APNEA DETECT

MODULE INFO

MODULE SETUP

ZERO


MODULE RESET

9.1 Introduction

ETCO₂(End-Tidal CO₂) is a device to see the concentration of end-tidal carbon dioxide, which uses a method of measurement based on the non-dispersed IR absorption of CO₂ using IR ray by sampling a certain part of respiration through pipe during respiration.

EtCO₂ connector position and accessory (Sidestream, Respironics)

EtCO₂ Connector



LoFlo sidestream CO₂ sensor and connector



Sidestream sensor



Sidestream sensor connector

BM3VET Touch Operation Manual

EtCO₂ accessories for sidestream applications


EtCO₂ monitoring accessory uses the accessories for LoFlo™ sidestream module of Respirationics Company.

The airway adapters for sidestream intubated applications			
3473ADU-00		Airway Adapter Kit w/ Dehumidification Tubing	Weight: 4.5 grams Deadspace – adds approximately 7 cc of deadspace Intended for use when monitoring Animals with ET Tube sizes >4.0 mm
3473INF-00		Airway Adapter Kit w/ Dehumidification Tubing	Weight: 5.8 grams Deadspace – adds approximately 1 cc of deadspace Intended for use when monitoring Animals with ET Tube sizes <=4.0 mm

BM3VET Touch Operation Manual

Connecting the LoFlo Sample Kit

1. The sample cell of the sampling kit must be inserted into the sample cell receptacle of the LoFlo CO₂ Module as shown in Figure 1. A “click” will be heard when the sample cell is properly inserted.




2. Inserting the sample cell into the receptacle automatically starts the sampling pump. Removal of the sample cell turns the sample pump off.

3. To remove the sampling kit sample cell from the sample cell receptacle, press down on the locking tab and pull the sample cell from the sample cell receptacle.

BM3VET Touch Operation Manual

EtCO2 connector position and accessory (Mainstream, Respironics)


EtCO2 Connector



CAPNOSTAT 5 mainstream CO2 sensor and connector



Mainstream sensor



Mainstream sensor connector

EtCO2 accessories for mainstream applications

EtCO2 monitoring accessory uses the accessories for CapnoStat 5 microstream sensor of Respironics Company.

BM3VET Touch Operation Manual

The airway adapters for mainstream intubated applications			
6063-00			Single-Animal Use Airway Adapter Intended for use when monitoring Animals with ET Tube sizes >4.0 mm
6312-00			Single-Animal Use Airway Adapter Intended for use when monitoring Animals with ET Tube sizes <=4.0 mm
7007-00			Reusable Airway Adapter
7053-00			Reusable Airway Adapter

BM3VET Touch Operation Manual

Connecting the CAPNOSTAT® 5 CO₂ Sensor to the Host System

1. Insert the CAPNOSTAT 5 CO₂ Sensor connector into the receptacle of the host monitor as shown in Figure 1.





Figure 1

2. Make sure the arrows on the connector are at the top of the connector and line up the two keys of the connector with the receptacle and insert.


3. To remove the connector, grasp the body portion of the connector back and remove.

Note: Do not remove by pulling cable.


Shown below is the CAPNOSTAT 5 CO₂ Sensor connection to a Respironics Novamatrix CO₂ adapter



Shown below is the CAPNOSTAT 5 CO₂ Sensor with an Animal circuit:



9.2 EtCO₂ Parameter Window



S: Display of apnea setting time in seconds

Upper/lower limit value of alarm: Display of alarm setting range value for concentration of CO₂

EtCO₂: Display of concentration value of carbon dioxide

AWRR: Display of the number of respirations per minute

FiCO₂: Display of concentration value of carbon dioxide during inspiration

Note

EtCO₂ waveform is always displayed if cable is connected.

9.3 EtCO2 Parameter Setting Menu

ALARM: A menu to set the alarm limit

WAVE SCALE: menu to set the size of waveforms on the on-screen

SWEEP SPEED: Speed to draw the signal waveform. (6.25mm/s, 12.5mm/s, 25mm/s)


APNEA DETECT: Menu for the detection of apnea

MODULE INFO.: Menu where you can see the MODULE information

MODULE SETUP: Menu to set module information.

ZERO: Atmospheric pressure and zero adjustment menu

MODULE RESET: EtCO2 MODULE menu to initialize the run




BM3VET Touch Operation Manual

ALARM LIMIT

(Upper/lower limit value of alarm)

Upper/lower limit value of alarm differs depending on the position of measurement.

The basic setting range of alarm setting value for EtCO₂, FiCO₂, AWRR, APNEA.



The following table shows standard alarm limit of parameter and setting value of scale when setting the label.

Parameter	Adult			Neonatal		
	Low	High	Scale	Low	High	Scale
EtCO ₂	0	98	40	0	98	40
FiCO ₂	0	20		0	20	
AWRR	0	100		0	100	
APNEA	0	40		0	40	

BM3VET Touch Operation Manual


WAVEFORM SCALE

(Measured waveform scale setting)

This sets the range of measured waveform versus pressure.

Selectable numerical value means the maximum pressure range value that is shown with waveform.


Pressing the knob switch key and then selecting the desired range value displays the selected pressure range value below the upper dotted line among two dotted lines in the left middle of wave window.



EtCO2 SWEEP SPEED

EtCO2 speed is 6.5mm/s.


Speed is changeable to 6.25, 12.5, 25mm/s.



BM3VET Touch Operation Manual


APNEA DETECT

Turn the APNEA detection alarm off and on



APNEA ALARM: This performs a function to set the display of apnea message alarm.

This displays a “apnea” message at the center of parameter window as shown in the figure below with apnea alarm on in case of apnea until the set apnea period is passed through.




With apnea alarm off, measured values are displayed instead of message.



BM3VET Touch Operation Manual

MODULE INFO

This is information for handling the EtCO2 module.



SENSOR PN(part number) : The sensor part number

OEM ID : The id is a 7bit identifier which is set at the factory to a unique value for each OEM.

SENSOR SN : The serial number of the module.

HW REVISION NUM : The hardware version number of the module.

TOTAL USE TIME : Total use time of the module.

LAST ZERO TIME : This is the total time that has elapsed with the sensor since the last zero.

PUMP TOTAL USE TIME : This is the total time the pump has been on.(LoFlo only)

PUMP MAX USE TIME : This value indicates the maximum rated lifetime of the sampling pump.
(LoFlo only)

BM3VET Touch Operation Manual

MODULE SETUP

This is information for handling the EtCO₂ module.

Setting	Value
BAROMETRIC PRESSURE:	760
GAS TEMPERATURE:	36.0
NO BREATH DETECT TIMEOUT:	20S
O2 COMPENSATION:	16
ANESTHETIC AGENT:	0.0
CURRENT ETCO2 TIME PERIOD:	20 S
CURRENT ETCO2 UNIT:	mmHg
BALANCE GAS:	ROOM AIR
SLEEP MODE:	NORMAL OP.
DISABLE SAMPLING PUMP:	NORMAL OP.
ZERO GAS TYPE:	ZERO ON ROOM AIR

BAROMETRIC PRESSURE:

This setting is used to set current Barometric Pressure.

GAS TEMPERATURE:

This setting is used to set temperature of the gas mixture. This setting is useful when bench testing using static gasses where the temperature is often room temperature or below.

NO BREATH DETECT TIMEOUT:

This setting is used to set the no breaths detected time-out. This time-out is the time period in seconds following the last detected breath at which the Capnostat will signal no breaths detected.

O2 COMPENSATION

ANESTHETIC AGENT

BALANCE GAS:

Use this setting to correct for the compensation of the gas mixture administered to the Animal. Anesthetic agent is ignored when the balance gas is set to helium.

BM3VET Touch Operation Manual

- CURRENT ETCO₂ TIME PERIOD:** This setting is used to set the calculation period of the ETCO₂ value. The end-tidal CO₂ value is the highest peak CO₂ value of all end of expirations (end of breaths) over the selected time period. If less than two breaths exist in the selected time period, the value will be the maximum ETCO₂ value for the last two breaths.
- CURRENT CO₂ UNIT:** Continuous waveform mode commands (the CO₂ Waveform Mode command [command 80h] and the CO₂/O₂ Waveform Mode command [command 90h]) MUST NOT be active when this command is used otherwise this command will be ignored and the setting will remain unchanged.
- SLEEP MODE:** Sleep mode is used to save power when the host monitor is in standby mode. There are two sleep modes available for the Capnostat. Using Sleep Mode 1 maintains the heaters so the Capnostat is able to run immediately after exiting the sleep mode. Mode 2 will require the Capnostat to go through its warm up sequence when exiting this mode and a delay will be introduced until the system has stabilized.
- ZERO GAS TYPE:** When performing a zero on room air, this setting should be set to room air (the default). Only change to nitrogen (N₂) when performing a zero on 100% N₂ gas; this is provided for use in a laboratory environment.
- DISABLE SAMPLING PUMP:** This setting allows the pump to be forced off. In Normal Operating Mode, the pump will be turned on when the sampling cell is connected and no pneumatic system errors are detected. In Pump Disabled Mode, the pump will remain off in all circumstances.


BM3VET Touch Operation Manual

ZERO

This function is used to initiate a Capnostat zero.

A zero is used to correct for differences in airway adapter types.

The Capnostat zero must be performed free of any CO₂.



1. Set the Host to the zeroing function.
2. Connect the CAPNOSTAT 5 CO₂ Sensor
3. Place the CAPNOSTAT 5 CO₂ Sensor onto a clean and dry CO₂ adapter that is exposed to room air and away from all sources of CO₂, including the ventilator, the Animal's breath and your own.
4. Start the adapter zero. The maximum time for a CAPNOSTAT zero is 40 seconds. The typical time for a zero is 15~20 seconds.


Note

For best result, connect the CAPNOSTAT 5 CO₂ Sensor to an adapter and wait 2 minutes before performing the Adapter Zero procedure.

BM3VET Touch Operation Manual

MODULE RESET

This performs a function to reset handling the EtCO2 module.



Warning

If defibrillation is performed while doing CO2 monitoring, remove the CO2 FilterLine from Animal Defibrillation without removing the FilterLine can result in serious electrical burn, shock, or injury due to electric discharge energy.

9.4 TROUBLESHOOTING

Following is a list of some of the messages that may appear on the monitor when monitoring CO₂. The message should clear when normal operating criteria are met or a solution is found.

* **SENSOR OVER TEMP**

- Cause : The sensor temperature is greater than 40°C
- Solution : Make sure sensor is not exposed to extreme heat(heat lamp,etc.)

* **SENSOR FAULTY**

- Cause: One of the following conditions exist : Capnostat Source Current Failure
EEPROM Checksum Faulty , Hardware Error
- Solution : Check that the sensor is properly plugged in. Reinsert or reset the sensor if necessary.

* **SENSOR WARM UP**

- Cause : Sensor under temperature , Temperature not stable, Source Current unstable
- Solution : This error condition is normal at startup. This error should clear when the warm up is complete.

* **CHECK SAMPLING LINE**

- Cause : This error occurs whenever the pneumatic pressure is outside the expected range.
- Solution : Check that the sampling line is not occluded or kinked. Replace the sample line

* **ZERO REQUIRED**

- Cause : Zero Required , Zero Error
- Solution : To clear, check airway adapter and clean if necessary. If this does not correct the error, perform an adapter zero. If you must adapter zero more than once, a possible hardware error may exist.

* **CO₂ OUT OF RANGE**

- Cause : The value being calculated is greater than the upper CO₂ limit(150mmHg)
- Solution : If error persists, perform a zero.

* **CHECK AIRWAY ADAPTER**

- Cause : Usually caused when the airway adapter is removed from the Capnostat or when there is an optical blockage on the windows of the airway adapter. May also be caused by failure to perform Capnostat zero when adapter type is changed.
- Solution : To clear, clean airway adapter if mucus or moisture is seen. If the adapter is clean, perform a Capnostat zero.

BM3VET Touch Operation Manual

Note

In the following monitoring conditions, the measured values may be inaccurate. Read the measured values carefully.

1. When using this in an environment of using nitrous oxide gas of high concentration
2. When using this in an environment where abrupt temperature change takes place
3. When using this in an environment with severely high humidity.

Caution

- The measured values may be inaccurate when using this equipment for Animals who have very fast or irregular respiration.
- When measuring CO₂ from the Animal under the anesthesia, check it when gas mixture comes in. Otherwise, the measured result values may be inaccurate.
- When using an anesthesia machine that uses a volatile anesthetic, CO₂ values may be inaccurate.

10. TEMPERATURE

10.1 Outline

10.2 Temperature Data Window

10.3 Temperature Data Setup


ALARM LIMIT
UNIT SELECT

10.1 Outline

This function is used to indicate the changes of resistance generated by the changes of temperature in numbers. The function involves the process of transferring the changes into electric signals.

Temperature Connector and Measuring Cable

Temperature Connector




Temperature Measuring Cable



Note

Temperature probe is correctly positioned and fixed to do not disconnect on the Animal. Temperature cable is attached to the monitor.

10.2 Temperature Data Window




Note

The minimum measuring time required to obtain accurate readings at the specific body site is at least 3 minutes.

10.3 Temperature Data Setup

ALARM : Temperature measurement alarm set


UNIT: Temperature measurement unit set



ALARM

Alarm menu provide ALARM LIMIT and ALARM.


Setting numeric value is 15.0°C ~ 45.0°C.



BM3VET Touch Operation Manual

UNIT SELECT

Able to select unit between °C, °F.



Warning

To measure the peripheral temperature, attach the probe to the ankle or palm.

If the patient sweats heavily or moves violently, fasten the pad with surgical tape.

NOTE

When the measuring site is exposed directly to air, the temperature may be lower than normal. It take about 20 to 30 minutes to reach the equilibrium temperature after attaching the sensor.

11. PRINT

11.1 Print

Printer and Heat Sensitivity Paper
Function and Setup Menu

11.2 Paper Change

11.1 Print

Printer and Heat Sensitivity Paper


A printer used to print data onto thermal paper.

Size of the thermal paper roll: 58mm wide x 38mm in diameter any thermal paper of same size can be used for the printer.


Side View of Printer



Function and Setup Menu



1. Press the PRINT Key for continuous printing.
2. Select Printing Speed 25, 50 mm/s.



3. Set up ALARM PRINT in the MORE menu to activate print on ALARM.



ALARM





ALARM PRINT




ON

BM3VET Touch Operation Manual

4. Data is printed in a selected wave form along with personal information of the Animal.
3 channels select 3 parameters to print.



BM3VET Touch Operation Manual



BM3VET Touch Operation Manual


PRINTER KEY

This menu is setup printing time delay in normal printing.

There are two menus for time configuration. One is Real-time, another is Delayed Time.

Real-time: This configuration makes printing out the newest data when the Printer Key is pushed.

Delayed time: This configuration makes printing out the data after 5 seconds from the Printer Key is pushed.




BM3VET Touch Operation Manual

PRINTER TIME

This is configuration of printed time in normal printing.

If the print out is not stopped in manual by PRINTER KEY, BM3VET TOUCH print out for setup time after starting print out with PRINTER KEY. The configuration of time could be setup with 4 types in CONTINUOUS, 10 sec, 20 sec and 30 sec. The configuration of PRINTER KEY(Real-time/Delayed time) is applied at print out with PRINTER TIME configuration.



If there is no print sheet, no paper icon of  appears.

BM3VET Touch Operation Manual

Thermal Paper Storage

To avoid fading of traces or deterioration, follow these precautions:

Note
These precautions apply to both unused paper as well as paper that has already been run through the printer.

- Store in cool, dark locations. Temperature must be below 27°C (80°F). Relative humidity must be between 40% and 65%.
- Avoid exposure to bright light or ultraviolet sources such as sunlight, fluorescent, and similar lighting which causes yellowing of paper and fading of tracings.
- AVOID CONTACT WITH: cleaning fluids and solvents such as alcohols, ketones, esters, ether, etc.
- DO NOT STORE THERMAL PAPER WITH ANY OF THE FOLLOWING:
 - carbon and carbonless forms.
 - non-thermal chart papers or any other products containing tributyl phosphate, dibutyl phthalate, or any other organic solvents. Many medical and industrial charts contain these chemicals.
 - document protectors, envelopes, and sheet separators containing polyvinyl chloride or other vinyl chlorides.
- DO NOT USE: mounting forms, pressure-sensitive tapes or labels containing solvent-based adhesives.

To assure MAXIMUM TRACE IMAGE LIFE, thermal paper should be stored separately in: manilla folders, polyester or polyimide protectors.

Plastic document protectors, envelopes, or sheet separators made of polystyrene, polypropylene, or polyethylene will not degrade thermal traces in themselves. However, these materials afford no protection against fading from external causes.

Paper manufacturers advise us that these thermal products should retain their traces when properly imaged and stored for about 3-5 years.

If your retention requirements exceed these guidelines, we recommend you consider alternate image storage techniques.

11.2 Paper Change


1

Open the window of the printer.



2

Insert the paper roll offered with the product into the printing unit. Place the roll in a proper way so that the printed paper can roll out upwards.



3

Press the printer window until it is properly shut. Incomplete closure may cause failure in printing.



12. MESSAGE LIST

Function	Message	Details
ECG	LEAD FAULT	Cable is not properly connected.
SpO ₂	CHEK PROBE LEAD FAULT	Probe is off animal. Cable is not properly connected.
RESP	LEAD FAULT APNEA	Cable is not properly connected. APNEA gives an alarm.
NIBP	INFLATION FAILURE CHECK CUFF OVER PRESSURE DEFLATION FAILURE CHECK CUFF OVER TIME CUFF PRESSURE MEASUREMENT ERROR PULSE TOO WEAK	Cuff hose is not properly connected. Cuff pressure is excessive. Cuff is bent, preventing deflation. Measure time exceeds the preset Level. Measure signal absent
EtCO ₂	MODULE OFF SENSOR WARMUP CHECK ADAPTOR CHECK LINE APNEA ZERO IN PROGRESS SENSOR FAULTY	Module is not properly connected. Sensor is initializing Adaptor is not properly connected. Tube is not properly connected. APNEA gives an alarm. Zeroing procedure when necessary. Sensor is not properly measured
TEMP	LEAD FAULT	Cable is not properly connected.
ALARM	ALARM VOL.OFF SILENCED ALARM PAUSE 5MIN	Alarm volume is off. Alarm key is pressed once Alarm key is pressed twice
TREND	NO ANIMAL DATA	No Animal's data input.
PRINT	No paper Icon	No paper in the printer
BATTERY	BATTERY LOW	The battery level is low, automatically power off within 5 minutes.

13. DEFAULT SETTING VALUE

13.1 HORSE-ICU Mode

Alarm level

	High	Medium	Low	Message
Asystole	0			
VTAC/VFIB	0			
VTAC	0			
PVC			0	
ST			0	
HR		0		
NIBP - S		0		
NIBP - M		0		
NIBP - D		0		
SpO ₂			0	
SpO ₂ -Rate				0
RR				0
RR-Apnea				0
T(° C)				0
EtCO ₂			0	
FiCO ₂				0
AWRR			0	
LEAD FAULT				0
LOW BATTERY				0

BM3VET Touch Operation Manual

Parameter Limits

	Low	High
HR	60	150
NIBP-S	80	200
NIBP-M	50	170
NIBP-D	30	150
SpO ₂	90	100
SpO ₂ -Rate	60	150
RR(RESPI)	15	100
RR-Apnea	0	20
T °C/° F	30.0/86.0	42.0/107.6
ST	-0.4	0.4
PVC	0	20
AWRR	10	30
EtCO ₂	25	50
FiCO ₂	0	5

BM3VET Touch Operation Manual

Display

Patient Age	30
Primary ECG	II
Arrhythmia	LETHAL
Detect Pace	Off
Print Waveform1	LEAD II
Print Waveform2	SpO2
Print Waveform3	Resp
Alarm Print	On
NIBP Interval	Off
NIBP Cuff Size	LARGE
RR(Resp) Lead	II
Alarm Volume	Off
QRS Volume	Off
Pulse Volume	Off
ECG Lead Fault	Message
SpO2 Check Probe	Message
Units for Height	cm
Units for Weight	kg
Temperature Units	° C
NIBP Limit Type	Systolic
ECG Filter	Monitor
PVC	ON
ST	ON

BM3VET Touch Operation Manual

13.2 DOG-ICU Mode

Alarm level

	High	Medium	Low	Message
Asystole	0			
VTAC/VFIB	0			
VTAC	0			
PVC			0	
ST			0	
HR		0		
NIBP - S		0		
NIBP - M		0		
NIBP - D		0		
SpO ₂			0	
SpO ₂ -Rate				0
RR				0
RR-Apnea				0
T(°C)				0
EtCO ₂			0	
FiCO ₂				0
AWRR			0	
LEAD FAULT				0
LOW BATTERY				0

Parameter Limits

	Low	High
HR	60	160
NIBP-S	80	200
NIBP-M	50	170
NIBP-D	30	150
SpO ₂	90	100
SpO ₂ -Rate	60	160
RR(RESPI)	15	100
RR-Apnea	0	20
T °C/°F	30.0/86.0	42.0/107.6
ST	-0.4	0.4
PVC	0	20
AWRR	10	30
EtCO ₂	25	50
FiCO ₂	0	5

BM3VET Touch Operation Manual

Display

Patient Age	30
Primary ECG	II
Arrhythmia	LETHAL
Detect Pace	Off
Print Waveform1	LEAD II
Print Waveform2	SpO2
Print Waveform3	Resp
Alarm Print	On
NIBP Interval	Off
NIBP Cuff Size	MEDI.
RR(RESP) Lead	II
Alarm Volume	Off
QRS Volume	Off
Pulse Volume	Off
ECG Lead Fault	Message
SpO2 Check Probe	Message
Units for Height	cm
Units for Weight	kg
Temperature Units	° C
NIBP Limit Type	Systolic
ECG Filter	Monitor
PVC	ON
ST	ON

BM3VET Touch Operation Manual

13.3 PUPPY-ICU Mode

Alarm level

	High	Medium	Low	Message
Asystole	0			
VTAC/VFIB	0			
VTAC	0			
PVC			0	
ST			0	
HR		0		
NIBP - S		0		
NIBP - M		0		
NIBP - D		0		
SpO ₂			0	
SpO ₂ -Rate				0
RR				0
RR-Apnea				0
T(° C)				0
EtCO ₂			0	
FiCO ₂				0
AWRR			0	
LEAD FAULT				0
LOW BATTERY				0

Parameter Limits

	Low	High
HR	70	180
NIBP-S	80	200
NIBP-M	50	170
NIBP-D	30	150
SpO ₂	90	100
SpO ₂ -Rate	70	180
RR(RESP)	15	100
RR-Apnea	0	20
T °C/° F	30.0/86.0	42.0/107.6
ST	-0.4	0.4
PVC	0	20
AWRR	10	30
EtCO ₂	25	50
FiCO ₂	0	5

BM3VET Touch Operation Manual

Display

Patient Age	30
Primary ECG	II
Arrhythmia	LETHAL
Detect Pace	Off
Print Waveform1	LEAD II
Print Waveform2	SpO2
Print Waveform3	Resp
Alarm Print	Off
NIBP Interval	Off
NIBP Cuff Size	SMALL
RR(RESPI) Lead	II
Alarm Volume	Off
QRS Volume	Off
Pulse Volume	Off
ECG Lead Fault	Message
SpO2 CHECK Probe	Message
Units for Height	cm
Units for Weight	kg
Temperature Units	° C
NIBP Limit Type	Systolic
ECG Filter	Monitor
PVC	ON
ST	ON

BM3VET Touch Operation Manual

13.4 CAT-ICU Mode

Alarm level

	High	Medium	Low	Message
Asystole	0			
VTAC/VFIB	0			
VTAC	0			
PVC			0	
ST			0	
HR		0		
NIBP - S		0		
NIBP - M		0		
NIBP - D		0		
SpO ₂			0	
SpO ₂ -Rate				0
RR				0
RR-Apnea				0
T(°C)				0
EtCO ₂			0	
FiCO ₂				0
AWRR			0	
LEAD FAULT				0
LOW BATTERY				0

Parameter Limits

	Low	High
HR	70	200
NIBP-S	80	200
NIBP-M	50	170
NIBP-D	30	150
SpO ₂	90	100
SpO ₂ -Rate	70	200
RR(RESP)	15	100
RR-Apnea	0	20
T °C/°F	30.0/86.0	42.0/107.6
ST	-0.4	0.4
PVC	0	20
AWRR	10	30
EtCO ₂	25	50
FiCO ₂	0	5

BM3VET Touch Operation Manual


Display

Patient Age	30
Primary ECG	II
Arrhythmia	LETHAL
Detect Pace	Off
Print Waveform1	LEAD II
Print Waveform2	SpO2
Print Waveform3	Resp
Alarm Print	Off
NIBP Interval	Off
NIBP Cuff Size	SMALL
RR(RESP) Lead	II
Alarm Volume	Off
QRS Volume	Off
Pulse Volume	Off
ECG Lead Fault	Message
SpO2 Probe Off	Message
Units for Height	cm
Units for Weight	kg
Temperature Units	° C
NIBP Limit Type	Systolic
ECG Filter	Monitor
PVC	ON
ST	ON

14. TROUBLE SHOOTING


14.1 Noise in ECG

- Gel is dry
- Electrodes does not stick well to skin




14.2 SpO2 malfunction


Connectors of the equipments are in bad condition?




14.3 Temp malfunction



14.4 NIBP malfunction




14.5 Abnormality in NIBP measurements




14.6 Failure in battery recharge


(the battery does not fully recharge in 6 hours or more)




14.7 Power failure



14.8 Periodic noises



14.9 Print failure



15. SPECIFICATION

Ease of use

Additional Function

Monitor Environmental Specifications

Power

Specification

Accessories Included

Option

BM3VET Touch Operation Manual

Ease of use

- Battery operation
- Attached printer
- Table and graphic trend

Additional Function

- LAN Connection

Monitor Environmental Specifications

- Operating Temperature : 15°C to 40°C (59°F to 104°F)
- Storage Temperature : - 10°C to 60°C (14°F to 140°F)
- Humidity : 20% to 95% RH
- Operating Attitude : 70(700) to 106Kpa(1060mbar)

Power

- AC 100-240V (50/60Hz)
- Adapter 18 V, 2.8 A

BM3VET Touch Operation Manual

Specification

Display, Resolution	8.0" color TFT, 800 x 600 pixels
Dimension, Weight	238(W) x 250(H) x 163(D) mm, Approx. 3.1kg
Parameter	ECG, Heart Rate, Respiration Rate, SpO2, Pulse Rate, Systolic BP, Diastolic BP, Mean BP, 1 x Temperature, EtCO2, FiCO2, Airway Respiration Rate
Trace	3 waveforms : 1*ECG, SpO2, RR or EtCO2 Sweep speed : 6.25, 12.5, 25, 50 mm/sec
Indicators	Categorized alarms (3 priority levels), QRS beep & SpO2 pulse beep, Percent(%) SpO2 pitch tone Battery status, External power LED, Touch screen, Rotary knob
Interfaces	DC input connector : 12 to 18VDC, 2.5A Defibrillator Sync. Output : - Signal Level : 0 to 5V pulse - Pulse width : 100 ± 10 ms LAN digital output for transferring data, Nurse call system connection - 0.3A at 125VAC – 1A at 24VDC DC output : 5VDC, 1A Max USB Barcode Scanner, USB & SD memory data storage
Battery	Rechargeable Li-ion battery, 1hours for continuous working
Thermal Printer	Speed : 25, 50mm/sec, Paper width : 58mm
Data Storage	128hours trends, 20cases of 10sec alarm waveform
Language	English, French, Spanish, Italian, Germany, Chinese, Russian, Czech, Bulgarian, Portuguese, Romanian, Hungarian, Turkish, Polish, Korean
ECG Performance	
Lead type	3-lead, 5-lead(option)
Lead Selection	3-lead : I, II, III 5-lead : I, II, III, aVR, aVL, aVF, V
ECG waveforms	3-lead : 1 channel 5-lead : 1 channels
Heart Rate Range	30 – 300 bpm
Heart Rate Accuracy	±1bpm or ±1%, whichever is greater
Sweep speed	6.25, 12.5, 25, 50 mm/sec
Filter	-Diagnosis : 0.05Hz - 150Hz -Monitoring : 0.5 – 40 Hz -Moderate: 0.5 – 25Hz -Maximum : 5 – 25 Hz
S-T segment detection range	-2.0 to 2.0 mV
Arrhythmia analysis	ASYSTOLE, VTACH, VFIB, ,PVC,
Pacemaker Detection Mode	Indicator on waveform display (user selectable)
Protection	Against electrosurgical interference and defibrillation
Respiration Performance	
Method	Thoracic impedance
Channel selection	RA-LA or RA-LL
Measurement range	5 – 120 Breath per minute
Accuracy	±1 Breath per minute

BM3VET Touch Operation Manual

Apnea alarm	Yes
SpO2 Performance	
Saturation range	0 to 100%
Saturation accuracy	70 to 100% ± 2 digits 0 to 69% unspecified
Pulse rate range	30 to 254 bpm
Pulse rate accuracy	± 2 bpm
NIBP Performance	
Method	Oscillometry with linear deflation
Operation Mode	Manual/Automatic/Continuous
Measurement range	Large Pressure : 20 to 260 mmHg Medium Pressure : 20 to 230 mmHg Small Pressure : 20 to 120 mmHg
Accuracy	Meets accuracy requirements of ANSI/AAMI SP10:1992 and 2002
Temperature Performance	
Measurement range	15 to 45 °C (59 to 113 °F)
Accuracy	± 1 °C
Compatibility	YSI Series 400 temperature probes
Sidestream CO2 (Option)	
Measurement range	0 to 150 mmHg, 0 to 19%
Accuracy	0-40mmHg ± 2 mmHg, 41-70mmHg $\pm 5\%$ of reading 71-100mmHg $\pm 8\%$ of reading, 101-150mmHg $\pm 10\%$ of reading
Respiration rate	2 to 150 breath per minute
Respiration accuracy	± 1 breath per minute
Mainstream CO2 (Option)	
Measurement range	0 to 150 mmHg, 0 to 19%
Accuracy	0-40mmHg ± 2 mmHg, 41-70mmHg $\pm 5\%$ of reading 71-100mmHg $\pm 8\%$ of reading, 101-150mmHg $\pm 10\%$ of reading
Respiration rate	0 to 150 breath per minute
Respiration accuracy	± 1 breath per minute

BM3VET Touch Operation Manual

Accessories Included

1. Main body of BM3VET Monitor	1 EA
2. 3-Lead Animal Cable (MECA3(AHA), MECE3(IEC))	1 EA
3. 3-Lead Animal Extension Cable	1 EA
4. NIBP extension hose (NBPCBL-400)	1 EA
5. Reusable small animal NIBP cuff (ICUFF-430)	1 EA
6. SpO ₂ extension cable (SPCBL-400)	1 EA
7. Reusable animal SpO ₂ probe (SPASENS-400)	1 EA
8. DC Power Adaptor with Power Cord (18VDC/2.8A, BPM050S18F02)	1 EA
9. Operator`s Manual	1 EA
10. Thermal roll Paper (PAPER-400)	2 Roll
11. Reusable Temperature Probe (Rectal/esophageal, TEMPESENS-430)	1 EA

Option

1. Sidestream EtCO ₂ Module (Respironics)	1 SET
2. Mainstream EtCO ₂ Module (Respironics)	1 SET
3. Sidestream EtCO ₂ airway adapter sampling kit	1 EA
4. Mainstream EtCO ₂ airway adapter	1 EA
5. 5-Lead Animal Cable(MECA5(AHA), MECE5(IEC))	1 EA
6. 5-Lead Animal Extension Cable	1 EA

Abbreviations and Symbols

Abbreviations and symbols which you may encounter while reading this manual or using the monitor are listed below with their meanings.

Abbreviations

		A
A	amps	
AC	alternating current	
ADT	adult	
ARRYTHM	arrhythmia	
ASYS	asystole	
Auto, AUTO	automatic	
AUX	Auxiliary	
aVF	left foot augmented lead	
aVL	left arm augmented lead	
aVR	right arm augmented lead	
		B
BPM	beats per minute	
		C
C	Celsius	
CAL	calibration	
cm, CM	centimeter	
		D
D	diastolic	
DC	direct current	
DEFIB, Defib	defibrillator	
DIA	diastolic	
		E
ECG	electrocardiograph	
EMC	electromagnetic compatibility	
EMI	electromagnetic interference	
ESU	electrosurgical cautery unit	
		F
F	Fahrenheit	
		G
g	gram	
		H
HR	heart rate, hour	
Hz	hertz	
		I
ICU	intensive care unit	
IBP	invasive blood pressure	
Inc	incorporated	

BM3VET Touch Operation Manual

		K
kg, KG	kilogram	
kPa	kilopascal	
		L
L	liter, left	
LA	left arm, left atrial	
LBS	pounds	
LCD	liquid crystal display	
LED	light emitting diode	
LL	left leg	
		M
M mean,	minute	
m	meter	
MIN,	min minute	
MM, mm	millimeters	
MM/S	millimeters per second	
MMHG, mmHg	millimeters of mercury	
mV	millivolt	
		N
NIBP	noninvasive blood pressure	
NEO, Neo	neonatal	
		O
OR	operating room	
		P
PED	pediatric	
PVC	premature ventricular complex	
		Q
QRS	interval of ventricular depolarization	
		R
RA	right arm, right atrial	
RESP	respiration	
RL	right leg	
RR	respiration rate	
		S
S	systolic	
sec	second	
SpO2	arterial oxygen saturation from pulse oximetry	
SYNC, Sync	synchronization	
SYS	systolic	
		T
Temp, TEMP	temperature	
		U

BM3VET Touch Operation Manual

V	precordial lead	V
V	volt	
V-Fib, VFIB	ventricular fibrillation	
VTAC	ventricular tachycardia	

W

X

X	multiplier when used with a number (2X)
---	---

Symbols

&	and
°	degree(s)
>	greater than
<	less than
-	minus
#	number
%	percent
±	plus or minus

PRODUCT WARRANTY

Product Name	Veterinary Monitor
Model Name	BM3VET TOUCH
Approval Number	
Approval Date	
Serial Number	
Warranty Period	4 years from date of purchase (2 years in Europe)
Date of Purchase	
Customer Section	Hospital Name : Address : Name : Phone :
Sales Agency	
Manufacturer	

* Thank you for purchasing BM3VET TOUCH

* The product is manufactured and passed through strict quality control and through inspection.

BM3VET Touch Operation Manual

International Sales & service

Bionet Co., Ltd. :

5F, Shinsegae I&C Digital Center 61 Digital-ro 31 gil,
Guro-gu, SEOUL 08375, REPUBLIC OF KOREA

Tel : +82-2-6300-6410 / Fax : +82-2-6499-7789 / e-mail: sales@ebionet.com

Website: www.ebionet.com

U.S.A sales & service representative

Bionet America, Inc. :

2691, Dow Ave, Suite B
Tustin, CA 92780 U.S.A.

Toll Free: 1-877-924-6638 / Fax: 1-714-734-1761 / e-mail: support@bionetus.com

Website: www.bionetUS.com

European sales & service representative

MGB Endoskopische Geräte GmbH Berlin :

Schwarzschildstraße 6
D-12489 Berlin, Germany

Tel: +49(0)-30-6392-7000 / Fax: +49(0)-30-6392-7011 / e-mail: sales@mgb-berlin.de

Website: www.mgb-berlin.de

BIONET CO., LTD.

Product Name: BM3 Vet Touch