

QUICK GUIDE AND SPECIFICATIONS

1 DUAL GAS MODULE

Dual Gas sidestream module is intended for monitoring of intubated and ventilated or non-ventilated small and large patients under anesthesia. The Dual Gas module is a multi-gas analyzer measuring carbon dioxide (CO₂) and one of five anesthetic agents with manual selection of the specific agent type. This is intended for use only with the BM3Vet Touch, BM7Vet, and BM Pro Series (BM3Vet Pro / BM5Vet Pro / BM7Vet Pro).




2 INSTALLATION OF THE WATER TRAP TANK




***** Uninstalling the water trap tank is the reverse order of the installation.**

QUICK GUIDE AND SPECIFICATIONS

3 INSTALLATION OF SAMPLE LINE AND AIRWAY ADAPTER



4 MOUNTING DUAL GAS MODULE



SLIDE THE MODULE ONTO THE BRACKET





QUICK GUIDE AND SPECIFICATIONS

5 CONNECTING THE MODULE TO THE DEVICE

The Dual Gas module will be automatically detected when it is connected to BM3Vet Touch, BM7Vet, or BM Pro Series (BM3Vet Pro / BM5Vet Pro / BM7Vet Pro).

****Please refer to the operation manual of BM3Vet Touch, BM7Vet, or BM Pro Series for details of the Dual Gas module operation.*



QUICK GUIDE AND SPECIFICATIONS

5

SPECIFICATIONS

CO ₂	Range:	0 – 100 mmHg; 0 – 13.3 kPa; 0-10% CO ₂ STPD (standard temperature and pressure dry)
	Accuracy:	± (0.2 vol% + 4% relative)
	Rise Time:	400 ms (average)
Anesthetic Agents	Gases:	Isoflurane, Enflurane, Halothane, Desflurane, Sevoflurane
	Range: Iso., Enf., Hal., Sevo	0 – 6%
	Range: Desflurane	0 – 18%
	Resolution	0 .01%
	Accuracy	± (0.15% vol% + 4% relative)
	Rise Time	450 ms (average)
Calibration		Factory calibrated
Power Up Time		30 sec
Delay Time		< 4sec
Dimensions	L x W x H	175mm x 50mm x 85mm
Power	Sensor only (not including sampling pump)	12V 0.15A
Respiration	Range: Accuracy:	0-100 breaths/min ±1 breath/min
Flow Rate		Range and accuracy 170 ml/min ±20ml/min (TBD)
Environment Condition	Operating Temperature: Storage Temperature: Ambient Humidity: Ambient Pressure:	15 – 35°C -5 to 50°C 15-95% RH 70kPa to 106kPa (525mmHg to 795mmHg)
Weight		Approx. 400g
Water Removal System		Water trap tank
Sampling Line		Anesthetic gas tolerant standard sample line
Application	Veterinary use only	Species: Canine, Feline, Equine