

bionet

WORLD'S FIRST

WIRELESS

WEARABLE VETERINARY ECG SENSOR



Bionet introduces the world's first wireless ECG vet monitor, giving you unmatched accuracy, flexibility, and improved workflow. **72 hour continuous, real-time monitoring** and data at your fingertips with the VEMO app, making VEMO the preferred and optimized solution in pre- and post-op care. **All wirelessly and instantly on your iOS and Android tablets and smartphones with the intuitive, easy to use VEMO app.**

vemo



Wearable & 100% wireless 72 hour continuous ECG monitoring sensor ideal for pre- and post-op care

Lightweight 7.5g | Compact 90.0 mm x 28.0 mm | Optimized for small to medium dogs & cats

Free iOS & Android app | Bluetooth connectivity

Intuitive, easy to use VEMO app lets you simultaneously monitor vital signs – ECG, HR, RR – of up to **6 pets**

Real-time data updates via WiFi anywhere in your clinic

EVERYTHING INCLUDED:



VEMO ECG Sensor



Charging Cradle



Charging Cable



Hydrogel Adhesive | 8 pcs



WIRELESS WEARABLE VETERINARY ECG SENSOR

SPECIFICATIONS

Size	28 mm × 90 mm × 7.9 mm / 7.5g
Parameter	ECG, Heart Rate (40 ~ 300BPM), Respiration (5 ~ 25 RESP)
Sampling Rate	128 Hz
Network	BLE Paired
Battery	72 hrs with Full Battery Charged
Splash, Water, and Dust Resistant	IP 25
OS	Android / iOS

VEMO TABLET & SMARTPHONE APP REQUIREMENTS




Requires iOS 11.0 or later




Requires Android SW revision 9+




VEMO APP DASHBOARD



BATTERY CHARGE INDICATOR







GREEN LED = FULL CHARGE



BLUE LED = CHARGING

INCLUDED ACCESSORIES

VEMO ECG Wireless Sensor	1ea	
Charging Cradle	1ea	
Charging Cable	1ea	
Hydrogel Adhesive	8ea	

SIMPLIFY YOUR WORKFLOW



6 Patient Simultaneous Monitoring



Real-time Data Updates




Bluetooth Connectivity



Simple Easy Pairing



Easy Data Sharing

CONTACT BIONET TODAY | (714) 734-1760 | sales@bionetus.com 

to see how VEMO wireless ECG monitoring can streamline your workflow

bionetus.com

bionet

Bionet America, Inc. 2691 Dow Ave. Suite B, Tustin CA 92780