



How to connect with Vemo

1. Install Vemo Link Windows Program on your PC. (Windows Only). Run Vemo app from your tablet or mobile phone. From the Vemo app, press the menu button (red circle shown in the picture) on top left corner.




2. You will see the Host IP address and Port number from your tablet or mobile phone.




3. Open the Vemo Link app from the PC. Make sure PC and tablet or phone are all connected with same network.



4. Click "Connect" and window will pop up. Enter the IP address and port number from the Vemo (tablet or mobile phone) and press "OK".



5. When connection is succeeded, it will cast whatever shown from the Vemo app to the PC Vemo Link. Click each patient window screen and you can save the report as "PDF".



6. You can press "+" and "-" to make it bigger or smaller.

