


BT-Link Next Software



All Brio XVet Series monitors include Bionet's revolutionary BT Link Next software, automating vital sign recording during surgeries and capturing abnormal events. Need more than 1 monitor for your hospital? BT Link Next allows you to network multiple Brio XVet monitors for seamless central monitoring from any Windows-based PC in the hospital. At the end of each case, the software will produce a PDF report of vital signs in custom intervals as well as waveform reports. Users can attach, import, or upload these reports to a hospital PMS system or digital anesthetic record.

BT-Link Next Mobile

Your on-the-go smartphone/tablet solution for seamless monitoring of Brio XVet devices from anywhere in your hospital. Compatible with iOS and Android devices, this app lets you keep a close eye on your monitors from anywhere in the hospital and notifies you during any alarm event. Stay connected and in control with BT-Link Next Mobile.



bionet

714-734-1760
Sales@BionetUS.com
BionetUS.com

bionet

Enhance Anesthetic Safety with the Power
of Brio XVet Multiparameter
Monitors from Bionet



A Symphony of Veterinary Monitoring Precision and Safety

- 3 Models to choose from (X3Vet, X5Vet, X7Vet)
- Capacitive touch screen, intuitive interface
- Automated vital sign recording and surgical event tracking
- Clinically validated, veterinary specific technology
- Sleek, lightweight design

This product is for animal use only



“Bionet makes great, reliable monitors with excellent customer service. We have 2 multi-specialty hospitals and intend to only purchase Bionet monitors from now on”

Care Center Vets - Cincinnati

Brio XVET Specs

Brio Series	Brio X3Vet	Brio X5Vet	Brio X7Vet
Screen	8" Color TFT LCD capacitive Touchscreen	12.1" Color TFT LCD capacitive Touchscreen	15.6" Color TFT LCD capacitive Touchscreen
Resolution	1024 x 600	1280 x 800	1366 x 768
Dimension, Weight	9.3 x 8.7 x 5.3 inches, approx. 4.4lbs	12.1 x 10.7 x 5.5 inches, approx. 6.8lbs	15.5 x 12 x 6.1 inches, approx. 9lbs
Standard Parameters	3 lead ECG, SPO2, RR, NIBP, TEMP	3 lead ECG, SPO2, RR, NIBP, TEMP	3 lead ECG, SPO2, RR, NIBP, TEMP
Optional Parameters	5/7 channel ECG, ETCO2, AG, 2 channel IBP	5/7 channel ECG, ETCO2, AG, 2/4 channel IBP	5/7 channel ECG, ETCO2, AG, 2/4 channel IBP
Product Warranty	4 Years	4 Years	4 Years

Brio Features/Benefits

- Durable design that can handle any clinical environment
- Suntech NIBP technology provides easy-to-obtain blood pressures with a high degree of motion tolerance.
- Includes 2 types of SPO2 probes for various patient and surgical case types.
- Get unrivaled portability with up to 2-3 hours of cordless monitoring.
- See clearly from a distance with the Brio X7Vet - the largest screen size in it's class (15.6").
- Included BT Link Next Software interfaces with the monitor to allow for automated vital sign recording.
- Attach vital sign and event PDF reports to a patient anesthetic record or PMS software for enhanced data management.

Capnography Options

Optimize patient safety with our advanced capnography technology. Our EtCO2 monitoring solutions offer seamless integration, providing accurate real-time readings during surgery and making it easy for clinical-decision making during procedures.



Capnostat Mainstream EtCO2 Sensor
Item Code: B-RMCO2
Warranty: 2 years



LoFlo Sidestream EtCO2 Sensor
Item Code: B-RSCO2
Warranty: 2 years



Dual Gas Sidestream EtCO2 Module
Item Code: DG-SENSOR
Warranty: 2 years