


Batteries & Transformers

COMPATIBLE DEVICES




BM Rechargeable
battery: LEAD-ACID
Item Code: BM-BAT-1

BM3: D1C0400001 ~ D1E08xxxxx(04.2002 ~ 08.2004)
BM3Vet: D1C0400001 ~ D1E08xxxxx (04.2002 ~ 08.2004)
BM3VetPlus: D2D0200001 ~ D2D07xxxxx(2003.02 ~ 07.2003),
BMP0308 ~ BMP0408(BMP0308 ~ BMP0408)
BM3Plus: D2D0200001 ~ D2D07xxxxx(2003.02 ~ 07.2003),
BMP0308 ~ BMP0408(BMP0308 ~ BMP0408)



Rechargeable battery:
Ni-MH (2line, each line is
5cell, 12V, 1800mA/h)]
Item Code: BM-BAT-2

BM3: D1E0900001 ~ D1F11xxxxx(09.2004 ~ 11.2005)
BM3Vet: D1E0900001 ~ D1F11xxxxx(09.2004 ~ 11.2005)
BM3VetPlus: BMP0409 ~ BMP0509, D2F1000001 ~ D2F11xxxxx
(09.2004 ~ 11.2005)
BM3Plus: BMP0409 ~ BMP0509, D2F1000001 ~ D2F11xxxxx
(09.2004 ~ 11.2005)




Rechargeable battery:
Ni-MH(10 cell, 1line,
12V, 2200mA/h)]
Item Code: BM-BAT-3

BM3: D1F1200001 ~ D1G03xxxxx(12.2005 ~ 03.2006), D1G04 ~
D1H0300001 - 57(04.2006 ~ 03.2007), D1H0300058 - 110 ~
D1I0300150(03.2007 ~ 03.2008)
BM3VetPlus: D2F1200001 ~ D2H02xxxxx(12.2005 ~ 02.2007)
BM3Vet: D1F12 ~ D1I03xxxxx(12.2005 ~ 03.2008)
BM3Plus: D2F1200001 ~ D2H02xxxxx(12.2005 ~ 02.2007)



BM Rechargeable battery:
Li-ION, 3cell, 1line,
10.8V, 2200mA/h
Item Code: BM-BAT-4

BM3: D1I0400001 ~UP TO DATE (04.2008 ~UP TO DATE)
BM3Vet: D1I0400001 ~UP TO DATE (04.2008 ~UP TO DATE)
BM3: DDL0700001 ~UP TO DATE(07.2011 ~UP TO DATE)
BM3Vet Next: DDL0700001 ~UP TO DATE(07.2011 ~UP TO DATE)
BM3VetPlus: D2H0300001 ~ D2I07XXXXX(03.2007 ~ 07.2008),
D2I1200001 ~D2M12xxxxx(12.2008 ~12.2012)
BM3Plus: D2H0300001 ~ D2I07xxxxx(03.2007 ~ 07.2008),
D2I1200001 ~ D2Mxxxxxx(12.2008 ~2012)
BM3Plus (Suntech): D3L0700001 ~ D3M12xxxxx(05.2010 ~
12.2012), D3N0200001 ~ UP TO DATE(02.2013 ~ UP TO DATE)
BM5Vet: D4G0900001 ~ D4I10xxxxx(09.2006 ~ 10.2008)
BM5Vet: D7J0400001 ~ D7K10xxxxx(04.2009 ~ 10.2010)
BM5Vet: D8K1100001 ~ UP TO DATE(11.2010 ~ UP TO DATE)
BM3Vet Touch: All units



Rechargeable battery:
Li-ION (EtCo2) 6 Cell,
2 Line, 11.1V, 4400mA/h
Item Code: BM-BAT-5

BM5Vet: D4G0900001 ~ D4I10xxxxx(09.2006 ~ 10.2008)



Rechargeable battery:
Li-Ion, 6Cell, 1Line,
10.8V, 4300mAh/46.44Wh
Item Code: BM-BAT-6

BM7Vet: DEP0400096 ~ UP TO DATE
BM7Vet Pro



Rechargeable Li-ion
battery, 3.6V,
4300mAh/15.48Wh
Item Code: OXW-BAT

Oxy9Wave Vet: O2P040090 ~ UP TO DATE



BM Rechargeable
battery: Li-ION 3 Cell,
1 Line, 11.1V, 2200mA/h
Item Code: BM1-BAT

BM1: DBJ1000001 ~ UP TO DATE(10.2009 ~ UP TO DATE)
BM1 (Suntech): DCL0800001 ~ UP TO DATE(08.2011 ~ UP TO DATE)



Optional Rechargeable
Battery:
Li-ion, 6Cell 2Lines, 10.8V,
4800mAh
Item Code: BM-BAT-7

BM7Vet: DEPO0400096 ~ UP TO DATE
BM7Vet Pro

COMPATIBLE DEVICES



INPUT: 100~240V, 50-60 Hz, 1.5A
OUTPUT: **+18 VDC**, 2.8A
(Doesn't include power cable (PC-A(V))
Item Code: TRANS

**BM3, BM3PLUS, BM3VET,
BM3VET TOUCH, BM5VET, BM7VET
BM PRO SERIES (BM3VET PRO
/ BM5VET PRO / BM7VET PRO)**



INPUT: 100~240V, 50-60 Hz, 1.0A
OUTPUT: **+15 VDC**, 2.0A
(Doesn't include power cable (PC-A(V))
Item Code: BM1 TRANS

BM1VET



INPUT: 100~240V, 50-60 Hz, 0.4A
OUTPUT: **+5 VDC**, 2.0A
(Doesn't include power cable (PC-A-OX))
Item Code: OXW-TRANS

OXY9WAVE VET



110V power cable
Item Code: PC-A(V)

TRANS, BM1VET TRANS



110V power cable
Item Code: PC-A-OXW

OXW-TRANS



INPUT: 110~240V, 47-63 Hz, 0.48A
OUTPUT: **+9 VDC**, 1.5A
(Doesn't include power cable (PC-A-OX))
Item Code: OXP-TRANS

Oxy9Vet Plus



BM Rechargeable
battery: Li-ION
3 Line, 10.8V, 2150mA/h
Item Code: BM-BAT-3&5

**BM3Vet Pro
BM5Vet Pro**