

VetPort-X IV

Veterinary Portable Dental X-ray System

 **Distributed by
Bionet America**



3.3 ^{lbs}
Light Weight


0.4 ^{mm}
Focal spot

70 ^{kV}
Tube Voltage

PortView
Digital Intra-Oral Sensor

Wireless

Transfer and Share X-ray Image Files



**Smart
Handheld X-ray
for smart veterinarians**

- **Compact Size**
- **Portable camera type**
- **Easy to control UI**
- **Superior image quality**

Compact

170 (W) x 138 (D) x 155 (H) mm

LCD

Full Color

Battery

Li Polymer
22.2VDC, 1,000mAh



VetPort-X IV

Portable Dental X-ray System

● Compact Design

ZEN-PX4 offers a convenient diagnostic option through its compact design and portability.

● LCD Display

View images instantly with a PORTVIEW Sensor!

● Wi-Fi

Send images via Wi-Fi to Genoray's software or App.

● Micro SD Card

Store images when on the go!

This Product is for Veterinary Use Only.



PortView

Digital Intra-Oral Sensor

● High Resolution

Automatic electronic gain control = compatible with your x-ray.

● Ergonomic

Vertical bitewings: robust and comfortable design.

● Reliable

Sealed cord: protected against fluids.

Technical Specifications & Dimensions

X-ray Output	Tube Voltage	70 kV
	Tube Current	2.5 mA
X-ray Tube	Anode Type	Stationary
	Anode Angle	12.5°
	Focal Spot	0.4 mm
X-ray Control	Exposure Time Set	0.05 - 1.6 sec
Battery Capacity	Type	Li Polymer
Weight	Main Body	3.3 lbs (1.5 kg)
Dimensions	Main Body	5.5" x 6.8" x 10" (14 x 17.3 x 25.4 cm)

* The specifications are subject to change without notice.

Technical Specifications & Dimensions

SPECIFICATIONS	SIZE 1	SIZE 2
Image Size (H x V)	20 x 30 mm	26 x 34 mm
Number of Total Pixels (H x V)	1,000 x 1,506 pixels	1,300 x 1,706 pixels
Pixel Size	20 x 20 µm	
Pixel Pitch (H x V)	20 µm	
Sensor Type	CMOS	
Dynamic Range	16 bit (0 to 5000)	
Interface	USB 2.0	

* The specifications are subject to change without notice.



Bionet America, Inc.

www.bionetus.com | 2691 Dow Ave. Ste.B, Tustin, CA 92780

877-924-6638

saels@bionetus.com