

bionet

WORLD'S FIRST

WIRELESS

WEARABLE VETERINARY ECG SENSOR



Bionet introduces the world's first wireless ECG vet monitor, giving you unmatched accuracy, flexibility, and improved workflow. **72 hour continuous, real-time monitoring** and data at your fingertips with the VEMO app, making VEMO the preferred and optimized solution in pre- and post-op care. **All wirelessly and instantly on your iOS tablets and smartphones with the intuitive, easy to use VEMO app.**

vemo



Wearable & 100% wireless 72 hour continuous ECG monitoring sensor ideal for pre- and post-op care

Lightweight 7.5g | Compact 90.0 mm x 28.0 mm | Optimized for small to medium dogs & cats

Free iOS app | Bluetooth connectivity

Intuitive, easy to use VEMO app lets you simultaneously monitor vital signs – ECG, HR, RR – of up to **6 pets**

Real-time data updates via **WiFi** anywhere in your clinic

EVERYTHING INCLUDED:



VEMO ECG Sensor



Charging Cradle



Charging Cable



Hydrogel Adhesive | 8 pcs

bionet
bionetus.com

CONTACT BIONET TODAY | (714) 734-1760 | sales@bionetus.com

to see how VEMO wireless ECG monitoring can streamline your workflow

vemo



WIRELESS WEARABLE VETERINARY ECG SENSOR

SPECIFICATIONS

Size	28 mm × 90 mm × 7.9 mm / 7.5g
Parameter	ECG, Heart Rate (40 ~ 300BPM), Respiration (5 ~ 25 RESP)
Sampling Rate	128 Hz
Network	BLE Paired
Battery	72 hrs with Full Battery Charged
Splash, Water, and Dust Resistant	IP 25
OS	iOS


VEMO TABLET & SMARTPHONE APP REQUIREMENTS




Requires iOS 11.0 or later



VEMO APP DASHBOARD



BATTERY CHARGE INDICATOR




GREEN LED = FULL CHARGE



BLUE LED = CHARGING

INCLUDED ACCESSORIES

VEMO ECG Wireless Sensor
Charging Cradle
Charging Cable
Hydrogel Adhesive



SIMPLIFY YOUR WORKFLOW



CONTACT BIONET TODAY | (714) 734-1760 | sales@bionetus.com

to see how VEMO wireless ECG monitoring can streamline your workflow

bionetus.com

Bionet America, Inc. 2691 Dow Ave. Suite B, Tustin CA 92780

bionet