


**BM7Vet  
ELITE**



**BM5Vet  
ELITE**



**BM3Vet  
ELITE**



# BMVet ELITE SERIES

MULTI-PARAMETER  
TOUCH SCREEN  
MONITORS

**UNCOMPROMISING PERFORMANCE.  
COMPREHENSIVE CONNECTIVITY.**

**The all-new, best-in-class Bionet BMVet Elite Series** line of multi-parameter touch screen monitors all feature advanced cardiac analysis – recognizing and alarming for the three most common lethal arrhythmias, and reliable and accurate SunTech non-invasive blood pressure technology. All BMVet Elite models are equipped with motion tolerant and low-perfusion resistant pulse oximetry measuring up to 450bpm. They're versatile and lightweight, making them ideal for dental and wet stations. Easily mount to a wall mount or rolling cart to save precious clinic space. All come standard with WiFi, built-in printer, and free BT-Link PRO and BT-Link Mobile data management software.

**And Bionet gives you the industry's  
best warranty – 4 years.**



**CONTACT BIONET TODAY  
(877) 924-6638**

**sales@BionetUS.com**

**See what the all-new BMVet Elite  
can do for your clinic.  
BionetUS.com**

# BMVet ELITE SERIES

## Multi-parameter Touch Screen Monitors



**UNCOMPROMISING PERFORMANCE. COMPREHENSIVE CONNECTIVITY.**

	<b>BM7Vet</b> ELITE 	<b>BM5Vet</b> ELITE 	<b>BM3Vet</b> ELITE 
<b>Standard Parameters</b>	ECG, SpO <sub>2</sub> , NIBP (with SunTech technology), Respiration Rate, Temperature	ECG, SpO <sub>2</sub> , NIBP (with SunTech technology), Respiration Rate, Temperature	ECG, SpO <sub>2</sub> , NIBP (with SunTech technology), Respiration Rate, Temperature
<b>Optional Parameters</b>	Mainstream EtCO <sub>2</sub> , Sidestream EtCO <sub>2</sub> , Dual Gas (Anesthetic Gas + EtCO <sub>2</sub> ), Invasive Blood Pressure (Dual Port), Advanced Arrhythmia Package (7 channel ECG w/ 5 lead cable)	Mainstream EtCO <sub>2</sub> , Sidestream EtCO <sub>2</sub> , Dual Gas (Anesthetic Gas + EtCO <sub>2</sub> ), Invasive Blood Pressure (Dual Port), Advanced Arrhythmia Package (7 channel ECG w/ 5 lead cable)	Mainstream EtCO <sub>2</sub> , Sidestream EtCO <sub>2</sub> , Dual Gas (Anesthetic Gas + EtCO <sub>2</sub> )
<b>Enhanced Connectivity</b>	Includes BT-Link PRO Includes BT-Link Mobile WiFi and wired LAN ready	Includes BT-Link PRO Includes BT-Link Mobile WiFi and wired LAN ready	Includes BT-Link PRO Includes BT-Link Mobile WiFi and wired LAN ready
<b>Screen</b>	12.1" Color TFT LCD Touchscreen	10.4" Color TFT LCD Touchscreen	8" Color TFT LCD Touchscreen
<b>Data Saving</b>	168 Hours Trend Data Stored	168 Hours Trend Data Stored	168 Hours Trend Data Stored
<b>Printer</b>	Standard, Built-In	Standard, Built-In	Standard, Built-In
<b>Data Management Software</b>	BT Link (Windows) BT Link Mobile (iOS, Android) USB Data Port	BT Link (Windows) BT Link Mobile (iOS, Android) USB Data Port	BT Link (Windows) BT Link Mobile (iOS, Android) USB Data Port
<b>Warranty</b>	4 Years Parts and Labor 	4 Years Parts and Labor 	4 Years Parts and Labor 

**CONTACT BIONET TODAY | (877) 924-6638 | [Sales@BionetUS.com](mailto:Sales@BionetUS.com)**  
**See what the all-new BMVet Elite can do for your clinic.**

Bionet America, Inc. 2691 Dow Ave. Suite B, Tustin CA 92780